



CE



## Operation and maintenance manual

Edition 1 from 2020-12-18

**YETI-** special wheelchair for children and teenagers



**ATTENTION!** Users with sight, reading or cognitive impairment are requested to report to the Manufacturer in order to obtain information for the user in an appropriate format.

### **Purpose of the Manual**

The operation and maintenance manual contains basic information necessary for the correct use of the Yeti stroller, including the principles of safe use, description of its construction, operation and adjustment, preparation of the stroller for operation, how to maintain and clean the stroller, and warranty conditions.

It is designed for persons who directly care of children with disabilities (assistants child - parents, guardians), which will be transported in a stroller. It is also intended for doctors and physiotherapists operating the stroller.

### **ATTENTION! Keep this manual for future reference!**

Keep the manual in an easily accessible place.

Manufacturer's contact details:

- Postal address: MyWam sp.j., ul. Szczecińska 10, 41-516 Chorzów, Poland
- contact phone number: +48 32 733 11 31
- e-mail address: [mywam@mywam.pl](mailto:mywam@mywam.pl)



**ATTENTION!** Before using the stroller, each User is obliged to read and follow this manual.

**ATTENTION!** This manual is an integral part of the product and must always be attached to the product being sold.

User and maintenance manual version: EN

**Edition: first of: 2020-12-18**

MyWamKupiec, Bartold, Angres sp.j. reserves the right to make technical and commercial changes to the content of the Operation and Maintenance Manual without prior notice. Each change will be determined by the date the instructions were last updated .

## Table of Contents

▪ Introduction .....	5
1. Intended use of the Yeti stroller .....	5
2. Basic information about the Yeti stroller .....	5
3. Quick start .....	6
▪ Safety .....	7
4. Important instructions and safety rules .....	7
5. Nameplate .....	12
6. Warning labels .....	12
7. Marking of symbols .....	12
▪ Product .....	13
8. General construction of the stroller .....	13
9. Configuration and use – stroller frame.....	15
9.1. Frame folding mechanism.....	15
9.2. Unfolding the stroller frame .....	15
9.3. Folding the stroller frame .....	16
9.4. Transport lock. Installing and removing the transport lock .....	16
9.5. Swivel front wheels, solid. Locking the direction of travel of the front wheels.....	16
9.6. Rear wheels, solid. Installation and removal of rear wheels on the axle.....	17
9.7. Rear wheel amortization system .....	17
9.8. Ergonomic handle for pushing the stroller. Mechanism and angle adjustment of the handle. ....	17
9.9. Central brake. Applying and deactivating the brake.....	18
9.10. Tilting foot.....	18
9.11. Frame adjustment in terms of seat depth .....	18
9.12. Adjusting the frame in terms of changing the seat angle .....	18
9.13. Frame adjustment in terms of backrest height.....	19
9.14. Adjusting the frame within the angle of the backrest .....	19
9.15. Footrest. Regulation and use. ....	20
10. Configuration and use - therapeutic seat with backrest.....	20
10.1. Stroller upholstery. Installation of upholstery on the stroller frame. Disassembly. ....	20
10.2. Seat upholstery adjustment.....	21
10.3. Backrest upholstery adjustment .....	21
10.4. Five-point belts .....	22
11. Additional equipment .....	22
11.1. Safety barrier .....	22
11.2. Therapeutic table .....	23
11.3. Abduction and stabilizing belts .....	23
11.4. Seat minimizing kit .....	24
11.5. The head restraint.....	24
11.6. Adjustable canopy with window .....	24
11.7. Travel bag.....	25
11.8. Shopping basket.....	25
11.9. Wedge.....	25

11.10.	Raincover .....	25
11.11.	Winter gloves .....	26
11.12.	Winter sleeping bag .....	26
11.13.	Mosquito net .....	26
11.14.	Sun umbrella .....	26
12.	Technical specification of the Yeti stroller .....	27
12.1.	User parameters .....	27
12.2.	Material composition of the stroller .....	27
13.	Delivery and tools .....	27
14.	Assembly and disassembly of the stroller .....	28
15.	Preparation for use .....	28
16.	Using the stroller .....	29
16.1.	Carrying the stroller .....	29
16.2.	Getting on and off the stroller .....	29
16.3.	Correct position in a stroller .....	30
16.4.	Riding in a stroller .....	30
16.5.	Moving on slopes .....	30
16.6.	Overcoming thresholds or curbs .....	30
16.7.	Transport in a car and other vehicles .....	30
16.8.	Cleaning and care .....	31
17.	Storage of the stroller .....	31
17.1.	Preparation of the stroller for long-term storage .....	31
17.2.	Corrosion protection .....	31
19.	Re-use .....	32
20.	Service and periodic inspections .....	32
21.	Disposal .....	33
22.	Compliance with requirements for medical devices .....	33
23.	Manufacturer .....	35
24.	Notes .....	35
25.	List of authorized service points .....	35
26.	Warranty Card .....	36
27.	Warranty conditions .....	36

## ■ Introduction

Thank you for purchasing special **Yeti** stroller for young people. We hope that it will be the cause of many wonderful moments with your child and will be helpful every day. We make every effort to ensure that our products meet your expectations and the highest standards. We encourage you to contact our employees and follow us on our website [www.mywam.eu](http://www.mywam.eu) and on social media. There you can find out many interesting things about our products.

### 1. Intended use of the Yeti stroller

- The Yeti Special Stroller is designed to transport a child with neurological, orthopedic or other conditions that hinder or prevent him from moving independently.
- The Yeti stroller is a medical device intended for transporting a child with paralysis or paralysis of a four-limb, three-limb, one-sided or crossed limb.
- The Yeti stroller has been designed to carry one child.
- The Yeti stroller, through the use of many regulations, allows comfortable and safe adaptation to the individual needs of the child.
- **The Yeti stroller is designed to carry a child whose weight does not exceed 75kg, height from 135cm to 175cm.** Statistically, they are children from 11 to 18 years old.
- When choosing a stroller, however, the most important are the individual physical and mental characteristics of the child, its dimensions and weight, housing conditions, environment and, above all, the individual's illness.
- Among the disease entities of a particular child, among others, there can be mentioned such diseases as: children cerebral palsy, meningitis, condition after spinal injury, condition after craniocerebral trauma, condition after stroke, paralysis of various origin and other dysfunctions damaging the locomotor and preventing independent movement of the child.
- Use the stroller as directed by a specialist.
- The Yeti stroller is intended for both indoor (e.g. apartment, rehabilitation center) and outdoor (e.g. park, pavement), on hard and level ground, on which all four wheels touch the ground and their contact with the ground is sufficient to properly drive the wheels.
- The Yeti stroller is intended to be controlled and propelled only by an assistant (parent, guardian).
- The user of the Yeti stroller is both the child who will be transported in it and the child's assistant (his parent, guardian, doctor, physiotherapist) who will transport the child in the stroller and look after him.
- Yeti stroller with proper care and maintenance maintains its quality for a period of 3 years and this is the expected period of its use.
- The risks associated with using the stroller are described in detail in the chapter: "Important instructions and safety rules".

### 2. Basic information about the Yeti stroller

The Yeti youth stroller has many options, e.g.:

- changing the angle of inclination of the backrest (see point 9.14)
- backrest height adjustment (see point 9.13)
- seat angle adjustment (see point 9.12)
- seat depth adjustment (see point 9.11 and 11.13)
- footrest adjustment (see point 9.15)
- head stabilization (see point 11.5)

- torso stabilization
- optional equipment

(see point 10.4, 11.3, 11.4, 11.9)

(see point 11)

### 3. Quick start










To quickly use the stroller it must be unpacked and prepared for use according to item 14 Assembly and disassembly of the stroller and item 15 Preparation for use. The following is a pictorial diagram, but it will not replace careful reading of the complete manual.



## ▪ Safety

### 4. Important instructions and safety rules

#### General tips:

-  • **WARNING!** Risk of accident - before using the stroller, the user is obliged to read and follow the entire contents of this manual.
-  • **WARNING!** Risk of accident - Before starting and while using the stroller, it is essential to comply with the safety instructions and rules described in this manual, and comply with the other instructions described in this manual.
-  • **WARNING!** Failure to follow the manual, in particular those regarding the moving parts of the stroller, may be dangerous for the user.
  - Always read the operating and maintenance manual before using the stroller.
-  • **WARNING!** Never leave your child unattended in a stroller, as there is a risk of injury and there is a risk of accident!
-  • **WARNING!** The product contains small parts that can be dangerous if swallowed by a child!
-  • **WARNING!** Risk of injury - always use all available security systems for the child using a stroller. These include: five-point belts, seat minimizing kit, frame wrappers, abduction and stabilizing belts, wedge, safety barrier or in its place a therapeutic table, footrest and straps for placing feet on the footrest.
-  • **WARNING!** There is a risk of burns when using the stroller - be careful when transporting a child in a stroller at high and low temperatures and in contact with the skin, e.g. in the sun, in cold weather, in a sauna, because the temperature of the stroller surface may rise when the stroller is exposed to solar radiation or the stroller surface temperature may drop when exposed stroller for frost.
-  • **WARNING!** When using the stroller or operating the stroller or folding the stroller, there may be a risk to the health of users and people in the immediate vicinity, caused by the dimensions of fixed and movable gaps and openings. Be particularly careful when doing this.
-  • **WARNING!** When using and operating the stroller, as well as when folding and adjusting the mechanisms, there may be a risk of entrapment and squeezing of parts of the body of the user or accompanying person in the openings and gaps between elements.
  - Be particularly careful when doing these operations.
  - After completing these operations, stabilize the positions of the individual components of the stroller by tightening the nuts, knobs and screws thoroughly. Each time you perform the above steps, you must first pull the child out of the stroller!
  - ATTENTION! Never exceed the maximum load.
  - ATTENTION! For Yeti stroller with load the maximum total load amounts to 80 kg. This value should not be exceeded.
  - ATTENTION! The maximum load of the travel bag is 2 kg.
  - ATTENTION! The maximum load of the shopping basket is 3 kg.
  - ATTENTION! The maximum therapeutic table load is 3 kg.
  - ATTENTION! The back pocket in the adjustable booth with the window is only used to place raincover or a mosquito net.
  - ATTENTION! Depending on the adjustments of the Yeti stroller - it concerns the position of the push handle, backrest height adjustment, backrest angle adjustment, seat depth adjustment, seat angle adjustment, footrest angle adjustment and footrest settings - this stroller may exceed the recommended maximum dimensions.






- **ATTENTION!** In the event of any medical incident related to the presence of a Yeti stroller, you must immediately notify the Manufacturer in writing by contacting via registered mail or email with confirmation of receipt using the addresses below:

Postal address: MyWam sp.j., ul. Szczecińska 10, 41-516 Chorzów, Poland






E-mail address: mywam@mywam.pl

- The therapy seat is not suitable for children under 10 years old
- It is forbidden to use the stroller as a ladder or as a base for climbing onto it to reach higher.
- The stroller must not be used in a way for which it has not been designed.

Tips in particular on construction and regulation:

-  • **WARNING!** All adjustments of the stroller and its elements, all settings, changes, fastenings and other activities related to matching the stroller should be carried out **WITHOUT A CHILD** placed in the stroller. Failure to do so may cause serious injury to a child and other people!
-  • **WARNING!** Adjustments made in the stroller, made incorrectly or not in accordance with the instructions in this manual, pose a serious danger to the child and other persons!
-  • **WARNING!** Due to the large adjustment possibilities of the stroller, when adjusting the seat angle and angle adjustment of the backrest, it is possible to set the position in which the child's head will be below the pelvis line. It is forbidden to use this position, as this could lead to the risk of aspiration or choking of the child.
- **ATTENTION!** Any change in the seat tilt changes the overall center of gravity, which affects the stability of the stroller.
- **ATTENTION!** Every change in the backrest's tilt causes a change in the overall center of gravity, which affects the stability of the stroller. Any change in backrest tilt requires five-point belt re-adjustment!
- **ATTENTION!** The central brake is not used to slow down the stroller during movement and should only be used to immobilize the stroller.
- The front steering wheels, if the direction of travel lock is off, may vibrate and lead to a sudden stop of the stroller and its falling forward when the surface is uneven. Therefore, with an uneven surface, it is recommended to block the front direction wheel lock and push the truck at a reduced speed.
- All adjustments and settings to be performed under the supervision and care of a doctor and a physiotherapist.

Information particularly regarding assembly and disassembly:

-  • **WARNING!** To avoid injury, make sure the child is at a far and safe distance when the stroller is assembled and disassembled.
-  • **WARNING!** The stroller must be fully assembled before use. In case of any problems at this stage, please contact the Seller and the Manufacturer urgently.
-  • **WARNING!** When mounting the stroller and when mounting the seat to the stroller frame, pay special attention to mount the seat in the designated place of the frame and secure all belts.
- Fastening the seat to other elements of the stroller frame is incorrect and may cause an accident!
-  • **WARNING!** It is forbidden to place additional weight on the stroller, which in particular applies to placing additional weight on the back of the back e.g. backpack or bag. This may cause the Yeti stroller to tip over and cause user injury! Does not apply to original additional accessories dedicated to the Yeti stroller not exceeding the permissible load.
-  • **WARNING!** When assembling and disassembling the stroller, keep a safe distance from moving parts of the stroller, as there is a risk of fingers, hands or other body parts being trapped.





- **WARNING!** When assembling and disassembling the wheels to the stroller frame, make sure and check that they are correctly fitted, as there is a risk of injury.
- **ATTENTION!** Do not unfold the stroller frame if the transport lock is in place, as it may damage it.
- All additional elements attached to the stroller affect its stability.

Tips in particular regarding preparation for use:



- **WARNING!** Before using the stroller, make sure that the stroller is properly assembled and all mechanisms, catches and fasteners are properly assembled, closed and prepared for use. In particular, this applies to: the mechanism for folding the frame, the mechanism for attaching the seat to the frame of the stroller and the installation of all four wheels.



- **WARNING!** Before riding the Yeti stroller, make sure that the child is secured and correctly adjusted and fitted with five-point belts, as there is a risk of injury.



- **WARNING!** The frame folding mechanism is a very important element of the construction of the stroller, so each time before using the stroller, pay special attention to the absence of cracks and damage to the folding mechanism and its proper operation. In the event of any suspicion you must not use the stroller, as this may endanger the user, child's assistant and other people! In this situation, it is recommended to urgently contact the Seller and the Manufacturer for advice.



- **WARNING!** Correct use of the frame folding mechanism is important for the safety of the child and other people.
- **ATTENTION!** The upholstery must always be dry before using the stroller.
- **ATTENTION!** Before each use of the stroller, the activities listed below must be carried out, and in the event of irregularities the stroller must not be used. In this situation, it is recommended to urgently contact the Seller and the Manufacturer for advice.

**The operations before each use of the stroller include:**



- Check the stroller for visible damage or other defects that could endanger the child during stroller operation.
- Check that the stroller frame has been unfolded correctly and the frame folding mechanism has been correctly closed (locked).
- Check that the seat has been correctly mounted to the stroller frame and that the seat attachment mechanism to the stroller frame has been properly closed (locked).
- Check the correct assembly of all four wheels.
- Check the good condition of the tread tires on all four wheels.
- Check good central brake operation.
- Correctly make all adjustments and adjust the stroller to the child's needs, including setting the appropriate depth of the therapeutic seat and height of the backrest, appropriate inclination of the therapeutic seat and inclination of the backrest, set the appropriate inclination of the footrest, adjust the appropriate length of five-point belts.
- Check that all adjustments are properly seated and secured.
- Check that the stroller upholstery is dry.
- Secure the child with five-point belts to prevent the child from falling out and sliding out of the stroller.
- Secure the child with all available stabilizing elements, such as five-point belts, seat minimizing kit, pelvic support belts, abduction and stabilization belts, wedge, footrest, footrest straps, safety barrier or therapeutic table. This is to protect the child from falling out and sliding out of the stroller and to protect the child's feet from falling from the footrest.

- **ATTENTION!** The stroller must not be used if the technical condition of the stroller or a single element of the stroller is bad or if the user has any doubts about it, the stroller is damaged or there is only a suspicion that damage has occurred or the stroller is not working properly. In this situation, it is recommended to urgently contact the Seller and the Manufacturer for advice.




Information particularly regarding the handling of the stroller:

- **ATTENTION!** It is forbidden to carry or lift the stroller with the child inside.
- **ATTENTION!** Due to the weight of the stroller and the inability to detach the seat from the frame, lifting the stroller is discouraged carrying it "in hands". If necessary, remove all accessories mounted additionally to the stroller and the rear wheels. The stroller should be lifted and carried by 2 people. Alternatively, you can lock the front wheels and pull the stroller behind you instead of carrying it.
- **ATTENTION!** Do not lift the stroller by the footrest or safety barrier, therapy table or any other additional accessories.




Information in particular regarding getting in and out of the stroller and the correct position in the stroller:

-  • **WARNING!** Do not stand on the footrest, as there is a risk of the stroller tipping over, injury to the user and risk of damage to the product.
-  • **WARNING!** Make sure that when getting in or out of the stroller, the child and no other person stand on the footrest, as this may cause damage and injury to the user!
- **ATTENTION!** When placing and removing a child from the stroller, the central brake must be activated.

Tips in particular regarding stroller driving:

-  • **WARNING!** The stroller must not be released by itself or pushed away, as it may cause an accident.
-  • **WARNING!** This product is not suitable for running or skating.
-  • **WARNING!** Use of the product for sport purposes is prohibited.
- **ATTENTION!** Make sure that during use of the stroller, the child's limbs do not get between the wheel spokes and the adjustable elements of the stroller.
- **ATTENTION!** When passing through a door or similar narrowing limiting the direction of the stroller, make sure there is adequate space on both sides of the stroller, which prevents your hand or other body parts from being trapped, as well as preventing damage to the stroller.
- When traveling on public roads, follow the Highway Code.

Tips in particular on moving on slopes:

-  • **WARNING!** It is forbidden to leave a child in the stroller without supervision of an adult and responsible assistant (parent, guardian), which in particular concerns leaving the child in the stroller on a slope, as it may cause an accident and injury to the user!
-  • **WARNING!** When using ramps, do not exceed the maximum load of the ramps, as this may cause injury to users and people around you.
-  • **WARNING!** When using ramps, make sure that the ramp is suitable, that the wheels of the stroller are high enough to overcome the driveway, that the stroller frame does not touch the driveway, and that help should be provided to another person as there is a risk of the stroller falling over, damage to the product and injury to the user.
- **ATTENTION!** Always remember that the location of the stroller's center of gravity changes when overcoming unevenness or slopes. In this case, be extremely cautious and secure the stroller.

Tips in particular on overcoming thresholds or curbs:



- **WARNING!** It is forbidden to climb or slide the stroller down the stairs, both with the child and the stroller itself, because it may cause an accident and injury to the user!

Information particularly regarding the transport of the stroller in a car and other vehicles:



- **WARNING!** It is forbidden to transport a child in a stroller and at the same time in any other means of transport, e.g. in a car, on a train, in a tram, on a bus, on a plane, on a ship or other, because it may cause injury.



- **WARNING!** It is forbidden to use the stroller seat as a car seat, as this may cause injury.



- **WARNING!** Risk of accident and injury to the user and other people - when transporting the stroller in other means of transport, e.g. in a car, bus, train, make sure that the stroller is correctly fastened to prevent injury to passengers during collision or sudden braking.
- **ATTENTION!** When transporting the stroller in other means of transport, e.g. in a car, on a bus, on a train, the stroller must be transported without a child and the stroller must be well protected against any damage.
- **ATTENTION!** The manufacturer is not liable for damages resulting from inadequate product protection during transport.
- **ATTENTION!** When transporting the stroller in other means of transport, e.g. in a car, on a bus, on a train, do not use one seat belt to secure both the passenger and the stroller at the same time.
- Use the transport lock when transporting the stroller and storing it.

Tips especially on maintenance and service:

- **ATTENTION!** The Manufacturer is liable for defects or technical damage to the stroller only caused by proper use of the stroller. The responsibility for the proper functioning of the stroller rests with the Manufacturer only if you strictly follow all the Manufacturer's instructions.
- **ATTENTION!** The Manufacturer is responsible for the proper functioning of the stroller only if original parts and original additional accessories are used.
- Only original spare parts that are supplied or recommended by the Manufacturer may be used in the stroller.
- Repairs to the stroller carried out by persons other than the authorized Service (see Chapter 19.) or the Manufacturer will void the warranty. The Manufacturer is not responsible for personal injury or damage to the product that may result from repairs carried out by persons other than the authorized Service center (see Chapter 19.) or the Manufacturer.
- The stroller should be cleaned and looked after regularly. Do not use abrasive cleaning agents.
- The Manufacturer is responsible for the proper functioning of the stroller only if it is purchased in a specialized medical store.
- The Manufacturer is not responsible for the appropriate selection of accessories for the given user.

## 5. Nameplate

- The nameplate contains basic information about the product: name, model, serial number (identification), maximum load, date of manufacture, confirmation of product compliance for medical devices, warning about proceeding in accordance with the Operation and Maintenance Manual and manufacturer's data (see Fig. 5.1).
- The nameplate is located on the crossbeam at the rear of the stroller frame, on its left side, glued from above.




Fig. 5.1

## 6. Warning labels

- The warning label contains important warnings for the user (see Fig. 6.1).

The warning label is located on the rear/side part of the frame on both sides.




Fig. 6.1

## 7. Marking of symbols

	Product serial number		The manufacturer's name
	Permitted total user weight	<p>DESCRIPTION "<b>WARNING!</b>" and the symbol in this manual are used to enhance the reader's attention to the content marked with this sign. Non-compliance with content marked with this description may endanger the health and safety of the user and accompanying persons, and may also damage the product.</p>	
	Date of manufacture		
	Conformity mark according to the Medical Devices Directive No. 93/42 EEC Annex VII		
	Warning about proceeding in accordance with the product Operation and Maintenance Manual		

▪ **Product**

**8. General construction of the stroller**

The Yeti Special Stroller for children consists of two main parts: the frame of the stroller with seat and backrest and very rich additional equipment that should be selected according to the individual needs of users.

• **Steel stroller frame with seat and backrest**

The basic elements of the stroller frame include:

- 1- Ergonomic handle for pushing the stroller - see point 9.8
- 2- String adjustment of excess upholstery - see point 10.3
- 3- Five-point belts - see point 10.4
- 4- Removable backrest and seat upholstery - see point 10.1
- 5- Slit for storing excess upholstery - see point 10.1
- 6- Frame wrapper - see point 10.1
- 7- Adjustable footrest with straps - see point 9.15
- 8- Tilting foot - see point 9.10
- 9- Central brake - see point 9.9
- 10- Rear wheel amortization system - see point 9.7
- 11- Rear wheels, solid - see point 9.6



- 12- Angle adjustment mechanism of the push handle - see point 9.8
- 13- Frame folding mechanism - see point 9.1
- 14- Backrest angle adjustment - see point 9.14
- 15- Seat depth adjustment - see point 9.11
- 16- Seat angle adjustment - see point 9.12
- 17- Backrest height adjustment - see point 9.13
- 18- Frame folding mechanism lock - see point 9.3
- 19- Angle adjustment mechanism of the entire footrest - see point 9.15
- 20- Backrest angle adjustment mechanism - see point 9.14
- 21- Transport lock - see point 9.4
- 22- Footrest length adjustment - see point 9.15
- 23- Quick rear wheels assembly mechanism - see point 9.6
- 24- Lock of the direction of travel of the front wheels - see point 9.5
- 25- Front steerable wheels, solid - see point 9.5

- **Optional equipment**

Its elements may include:

- Safety barrier - see point 11.1
- Therapy table - see point 11.2
- Abduction and stabilizing belts - see point 11.3
- Seat minimizing kit - see point 11.4
- Headrest with adjustment ⇕ - see point 11.5
- Adjustable canopy with window - see point 11.6
- Travel bag - see point 11.7
- Shopping basket - see point 11.8
- Wedge - see point 11.9
- Raincover - see point 11.10
- Winter gloves - see point 11.11
- Winter sleeping bag - see point 11.12
- Mosquito net - see point 11.13
- Sun umbrella - see point 11.14





## 9. Configuration and use – stroller frame

### 9.1. Frame folding mechanism

The frame folding mechanism is located inside the frame and on the frame in two places in its right part and in three places in its left part. Describing these places from above, from the side of the handle to push the stroller, these are:

- right and left bend - an element that allows folding the upper frame to reduce its dimensions
- mechanism of folding the frame with a lock



Fig. 9.1

### 9.2. Unfolding the stroller frame

To unfold (open) the stroller frame: lay the stroller frame horizontally on the ground and pull off the transport lock (Fig. 9.2.1). Then gently pull the backrest to the upright position (Fig. 9.2.2).



Fig. 9.2.1



Fig. 9.2.2

To set the stroller, pull the handle carefully by lifting it slightly up while bringing the stroller to the desired position (Fig. 9.2.3) until the frame folding lock engages (Fig. 9.2.4). The frame of the stroller will not unfold properly until the frame folding mechanism has been properly closed. This is evidenced by the clear "click" sound coming from the well-closed mechanism of folding the frame.



Fig. 9.2.3



Fig. 9.2.4



### 9.3. Folding the stroller frame

To fold (close) the stroller frame: release the frame folding mechanism lock by pulling it up and pull both frame folding handles up until the frame is released. (Fig. 9.3.1) Then fold the stroller frame folding it forward in the direction of the front wheels. Make sure that the push handle and footrest are in a non-colliding position. Then pull the backrest release mechanism. Put on a transport lock.



Fig. 9.3.1

### 9.4. Transport lock. Installing and removing the transport lock

The stroller is equipped with a transport lock, which serves to protect the folded (closed) stroller frame against its automatic opening during transport and storage.

The transport lock consists of two parts: an eyelet (female element) and a pin (male element). In order to secure the frame against its automatic opening, first fold the frame, and then put the locking eye on the locking pin so that both of these elements closely overlap (Fig. 4.4.).

To unlock the transport lock, before unfolding the frame, remove the locking eye from the locking bolt and only then can the unfolding of the stroller frame be started.



Fig. 9.4

### 9.5. Swivel front wheels, solid. Locking the direction of travel of the front wheels

The used, durable wheels guarantee the highest quality of stroller driving, driving pleasure and elegant appearance. The tires are tubeless, filled with a special gel that has shock absorbing properties very similar to those of inflatable wheels. The wheels are factory-fitted to the front wheel forks, which are finished with top steering axles from above.

Both front wheels can be set for straight-ahead or unlocked for turning, which makes it much easier to maneuver the stroller. Thanks to the locking of the direction of travel of the front wheels and the possibility of putting it on or taking off, there is clearly greater comfort for both the child and the assistant when driving the stroller.

It is recommended to remove the lock and allow the front wheels to turn while driving the stroller on a completely smooth surface, e.g. in shopping centers, school.



In the event of uneven surfaces, moving on sloping surfaces, when crossing thresholds or curbs, it is recommended to install the lock and block the turning of the front wheels.

To prevent the front wheels from turning, put the lock knob into the gap - see Fig. 9.5.1. To unblock the turning of the front wheels, pull back and turn the locking knob and leave them outside the gap - see Fig. 9.5.2.



Fig. 9.5.2

### 9.6. Rear wheels, solid. Installation and removal of rear wheels on the axle

The used, durable wheels guarantee the highest quality of stroller driving, driving pleasure and elegant appearance. The tires are tubeless, filled with a special gel that has shock absorbing properties very similar to those of inflatable wheels.

Thanks to the possibility of dismantling both the wheels and their axle, the trolley is smaller after folding, which facilitates its transport and storage.

The wheel assembly mechanism is mounted inside the rear wheels and ensures quick assembly and disassembly of the rear wheels. It is started by pressing the quick assembly clamp (Fig. 9.6.1) towards the center of the wheel, after which the rear wheel can be mounted on the rear wheel axle or the rear wheel can be removed from the rear wheel axle.

- Wheels mounting on axle

When installing the rear wheel on the rear wheel axle, press the clamp in the center of the wheel and put the wheel on the beginning of its axle. Then stop pressing the clamp and continue sliding the wheel onto its axle until it clicks. The correct wheel placement is demonstrated by the clear "click" sound coming from the quick axle rear axle mounting mechanism.

- Dismantling the wheels from the axle

In order to remove the rear wheel from the rear wheel axle, press the clamp in the middle of the wheel and slide the wheel out of its axle. - Fig. 9.6.2.



Fig. 9.6.1



Fig. 9.6.2

### 9.7. Rear wheel amortization system

The stroller is equipped with spring-loaded cushioning of the frame and rear wheels, located at the ends of the rear legs of the frame, when in contact with the crossbeam. Amortization works automatically. The adjustment is done by activating the up-down lever located from the bottom of the shock absorber. In order to facilitate the switchover of settings by the lever, simultaneous pressure on the shock absorber spring is allowed. When the lever is shifted, the shock absorber spring will be compressed (option for hard ground) or extended (soft ground).




Fig. 9.7

### 9.8. Ergonomic handle for pushing the stroller. Mechanism and angle adjustment of the handle.

The handle is used to push the stroller, providing a firm grip in various driving conditions. It has the ability to adjust the height of the handle for the child's assistant who drives the stroller - Fig. 9.8

- Adjustment mechanism.

These are two buttons on the outside of the frame, located on the joints of the handle to push the stroller.

- Angle adjustment of the stroller push handle

To change the angle of the stroller push handle, push both buttons of the angle adjustment mechanism of the push handle and change its angle. After setting the handle in the correct position, stop pressing the buttons, which will automatically lock the push handle in the desired position.



Fig. 9.8

### 9.9. Central brake. Applying and deactivating the brake



Fig. 9.9

The central brake pedal (Fig. 9.9), which when pressed blocks the two rear wheels of the stroller directly, is located in the middle of the crossbeam, at the rear of the stroller frame.

- Central brake on and off

To apply the central brake, press the brake pedal down towards the ground until it clicks. The correct activation of the brake is indicated by a clear "click" sound from the central brake.

To disengage the central brake, pull the brake pedal up.

### 9.10. Tilting foot

The tilting foot (Fig. 9.10) is located on the right side of the crossbeam, at the rear of the stroller frame. In order to use it, press the tilting foot and at the same time pull the handle to push the stroller, which in turn will cause the front of the stroller to rise up together with the front wheels.

Using the tilting foot makes it easier for the stroller to climb onto small steps or curbs.



Fig. 9.10

### 9.11. Frame adjustment in terms of seat depth



Fig. 9.11

The stroller has the option of adjusting the seat depth. To do this, activate the central brake so that the stroller is stable and immobilized, then loosen the screw with an Allen key from both sides of the frame (locking position) and extend or slide the seat frame (seat frame) evenly into the desired position. The correct position must be secured by tightening the screw on the sides of the frame. The frame can be pulled out leaving a minimum of 15cm overlaps of steel frame elements. See Fig. 9.11

The frame adjustment is in the range of 30 to 45 cm. The adjustment range of seat upholstery is described in point 10.2.

### 9.12. Adjusting the frame in terms of changing the seat angle.

The stroller has the ability to adjust the seat angle in the range of 10° to 30° (from level), every 5°. The seat angle adjustment mechanism consists of a comb on the right and left side of the stroller, to which the stroller frame is attached with a screw.

To set the right seat angle, activate the central brake so that the stroller is stable and immobilized, then unscrew the screw on both sides of the comb, set the seat according to the selected position and tighten the screw in the correct position. A Phillips screwdriver is required for this operation. See Fig. 9.12.

The upholstery adjustment range is described in point 10.3



Fig. 9.12

### 9.13. Frame adjustment in terms of backrest height

The stroller has adjustable backrest height. To do this, activate the central brake so that the stroller is stable and immobilized, then loosen the screw with the Allen key from the sides of the frame (locking position) and extend or slide the backrest (backrest frame) evenly into the desired position. Then secure the selected position by tightening the screw. The frame can be pulled out leaving a minimum of 15cm overlaps of steel frame elements. The frame adjustment is in the range of 65 to 85 cm. See Fig. 9.13

The upholstery adjustment range is described in point 10.3.



Fig. 9.14

### 9.14. Adjusting the frame within the angle of the backrest

The stroller has the ability to adjust the angle of the backrest, i.e. the backrest tilt from a sitting position to a lying position. This is done by means of a mechanism that consists of: a release cable, a left and right comb (element with an 8-point scale), and pins that fix the position on the comb connected by a cable. (See Fig. 9.14.1).



Fig. 9.14.1




Fig. 9.14.2

To change the backrest angle, activate the central brake so that the stroller is stable and immobilized, then pull the release cable up (see Fig. 9.14.2) and adjust the backrest to the desired value. The backrest should be locked by releasing the cable releasing the lock and making sure that the lock clicks on both sides of the stroller at the same time (it goes deep into the combs).

When changing the position of the backrest, it is necessary to adjust the upholstery. The upholstery adjustment range is described in point 10.3.

The angle adjustment mechanism of the backrest changes the position of the backrest on the back relative to the therapeutic seat in six positions, from sitting at right angles to almost fully lying position. The angle between the backrest and the seat is in the range of 90°-160° every 10° (measured from level).

Caution! Each change in the angle of the backrest requires a readjustment of the five-point belt length.



## 9.15. Footrest. Regulation and use.

The adjustable footrest is included in the standard equipment of the Yeti stroller. It consists of a footrest plate bolted to the footrest frame, with the possibility of moving the plate up or down in 7 positions to adjust the height of the footrest to the length of the child's thighs, in the range from 28.5cm to 43.5cm.



Fig. 9.15.1

- Footrest height adjustment.

To change the height of the footrest, unscrew the knobs on the two screws that secure the footrest to the frame on both sides of it, then set the footrest at the appropriate height and screw it back with both screws, tightening them with the knobs until you feel resistance. (Fig. 9.15.1)

- Angle adjustment mechanism of the entire footrest

The mechanism is located in the front of the seat and connects them to the footrest frame. They are two plastic, black joints with buttons inside (see Fig. 9.15.2). This

mechanism is used to change the angle of the footrest in relation to the seat, allowing you to choose the appropriate one of eight positions (range from 10° to -90°, every 20°, measured from level). To change the angle of the footrest, press the footrest angle adjustment mechanism on both sides of the footrest, then set the appropriate position and release the buttons. The footrest angle adjustment mechanism during use should be secured with upholstery wrappers (see point 10.2).



Fig. 9.15.2

- Use of a footrest.

It is recommended to position the footrest at such a height that the child's feet can rest freely on it with their entire surface, while the child's thighs are parallel to the seat surface. There are straps in the footrest plate to secure the child's feet. It is recommended to put them on when using the stroller.

Normally the footrest plate is in a horizontal position. However, during transport, storage, rehabilitation, getting on and off the Yeti stroller, the footrest plate must be folded and turned towards the footrest frame.

## 10. Configuration and use - therapeutic seat with backrest

### 10.1. Stroller upholstery. Installation of upholstery on the stroller frame. Disassembly.

The stroller upholstery consists of several elements: upholstery base, seat upholstery, back upholstery and frame wrappers.

- Upholstery installation.

Before adjusting the upholstery, the frame: backrest height, seat depth and footrest height must be adjusted accordingly. See point 7.12.

In order to mount the upholstery on the frame of the stroller, place the **upholstery base** in the place of the seat with the help of adjustable tapes (Velcro), stretching them thoroughly. The top of the upholstery base is characterized by longitudinal Velcro, to which the right upholstery will be later attached. Mount the upholstery base strips from the bottom, wrapping them on the seat frame (from the outside). The first Velcro from the footrest side is still mounted on the plastic fragment of the bent,

but from the inside of the seat. We wrap the excess of upholstery in the minimum seat setting at the back of the seat. The upholstery base should be properly stretched, allowing the wrapper to be attached to the frame in the places of the seat folds (in contact with the footrest, on the right and left).

**Frame wrappers** are designed to protect the places of the folds from friction of the user. To ensure comfort for the user, they should be mounted with Velcro down on the plastic elements of the bents.

We put the **backrest upholstery** on the base of the upholstery and fix it with straps at the back, wrapping the frame from the outside. If the minimum backrest height is used, the excess upholstery is wound from the top onto the frame and pulled with straps adjusting to the right distance through the lower buckle. If the backrest is extended to the maximum, we put the upper part on the frame and pull it through the upper buckle.

To properly mount the side covers of the stroller upholstery, it is recommended that the backrest be positioned to position 45° and secured with Velcro on the right and left of the stroller in two places on the side frame (when the frame folds). After placing the backrest in an upright position, the excess upholstery from the sides can be collected with strings with clamps.

Finally, put on the **seat upholstery**, arranging and adjusting it to the Velcro on the upholstery base. Excess upholstery drag under the upholstery of the backrest and wrap up (on the back) under the backrest upholstery strips. The belt that installs the seat upholstery (mounted underneath the seat embracing the frame from the outside) should fit over the first Velcro of the base upholstery and should not interfere with the frame wrapper. Attach the bottom part of the seat upholstery with Velcro to the footrest frame.

The stroller seat upholstery has an additional strip of material connecting it to the barrier. If it is not needed, it can be hidden by dragging in a dedicated slot - see Fig. 10.1.



Fig. 10.1

- Removing the stroller upholstery

We remove the upholstery using the previous steps in reverse order, i.e. unfasten the seat upholstery straps and remove this upholstery by "pulling" it away from the Velcro. Detach the sides of the backrest in the backrest upholstery and remove the upper backrest protection. Unhook all the straps and remove the backrest upholstery. Then remove the wrappers from the frame. The last step is to open all upholstery base straps and remove them from the frame. The upholstery must be protected against damage.

### 10.2. Seat upholstery adjustment

The seat can be adjusted in the range of 30 to 45cm. This is done by extending the frame (see point 9.11) and adjusting the upholstery to the required parameters.

- To change the seat depth parameters, remove the seat upholstery, adjust the seat frame in terms of depth by sliding in or out the seat frame to the desired size. The seat frame can be extended by 15 cm. Then place the seat upholstery on the seat base and attach it with Velcro (see item 10.1).

### 10.3. Backrest upholstery adjustment

The backrest can be adjusted from 65 to 85cm. This is done by extending the frame (see point 9.13) and adjusting the upholstery to the required parameters.

- In order to change backrest height parameters, detach the upper part of the backrest upholstery from the clamping straps at the back of the stroller and then adjust the upper frame of the stroller by extending or sliding the upper part of the backrest to the correct position. The adjustment range is 20cm. After choosing the right position, arrange the upper part of the upholstery adapting it to the frame. If the stroller is in the minimum backrest position, the upholstery may need to be rolled

up on the upper part of the frame. The upholstery from above should be fastened with straps to the buckles (lower in the case of a stroller with the smallest backrest, upper buckles in the case of the seat in the largest position).

- In order to correctly position the sides of the upholstery, the strings should be adjusted with the "pull" upholstery. For a stroller a sitting position, excess upholstery should be removed with strings in two places on each side. However, before setting the backrest to the lying position, the strings should be completely loosened, freeing the possibility of maximum unfolding of the upholstery.

#### 10.4. Five-point belts



Fig. 10.4.1

The Yeti stroller is standard equipped with professional five-point seat belts (Fig. 10.4.1). The straps have already been fitted to the stroller at the factory and should not be removed, as they are used to transport the child safely, prevent it from slipping, and are also helpful with additional stabilization of the child's pelvis and torso.

The buckle used to fasten them (see Fig. 10.4.2) facilitates getting in and out of the child. Sliding regulators, which are responsible for adjusting the length of the belts, are located on the right and left shoulder straps and on the fastening of the buckle. The belts are additionally equipped with a buckle fastening at chest height. This is an important element that helps to keep the belts in the desired position.

- Five-point belts adjustment. Adjust the five-point belts to your child before attaching them. To adjust the length of the belts, move the adjusters on the braces up or down and move the protectors accordingly. To adjust the height of the braces to be attached to the backrest, move the braces from low to high, or vice versa.

- Use. In order to put on the child five-point belts, put the braces on the right and left of the child and then attach both braces to the belt buckle, which is fixed on the seat. In order to remove five-point belts from a child, press the button on the upper (chest) and lower buckles and remove both braces from the buckle and slowly remove these braces from the right and left side of the child.



Fig. 10.4.2

Attention! Each change in the angle of the backrest requires a readjustment of the five-point belt length.

### 11. Additional equipment

#### 11.1. Safety barrier

The Yeti stroller can be equipped with a safety barrier that prevents the sitting child from falling forward in certain indications.

The safety barrier is mounted with the same fixings as the therapeutic table, so you cannot use both of these elements at the same time. Fastening the therapeutic table or safety barrier to the frame are two plastic elements (latches) mounted in the place of narrowing of the frame (see Fig. 11.1).

To install the safety barrier, open the two clamps on the railing and slide them onto the frame at the same time, latch. The barrier will be properly mounted only when the barrier covers the frame and snaps into place. This is evidenced by the clear "click" sound from a well-connected barrier attachment.



Fig. 11.1



To disassemble the barrier pull both barrier latches on the frame and pull out the barrier at the same time.

### 11.2. Therapeutic table

The Yeti stroller can be equipped with a non-slip therapeutic table with a maximum load of 3kg. The therapeutic table is mounted with the same fixings as the safety barrier so you cannot use both of these elements at the same time.

The attachment of the therapeutic table or safety barrier to the frame are two plastic elements (latches) mounted in the place of narrowing of the frame (see Fig. 11.2).

To assemble the table, open two clamps (latches) and simultaneously slide them onto the frame, latch. The table will not be properly installed until the latch engages the frame and snaps into place. This is evidenced by the clear "click" sound coming from a well-connected therapeutic table mount.

To disassemble the table, pull both latches on the frame and pull out the table at the same time.



Fig. 11.2

### 11.3. Abduction and stabilizing belts



Fig. 11.3.1

Abduction and stabilizing belts (Fig. 11.3.1) are used for additional pelvic stability and prevent the child from slipping out of the stroller. The correct way to use them is to fasten the child's thighs so that his pelvis adheres as close as possible to the backrest. This will ensure a healthier and more comfortable position of the child in the stroller.

The Abduction and stabilizing belts are mounted to the therapeutic seat in four places: two rear belt straps, attached to the rear seat frame with clamps (See Fig. 11.3.2), while the two front belt straps are attached to the trolley frame under the rail narrowing (sideways) by passing the belt through the factory loop in the tape (See Fig. 11.3.3)

- Adjustment. To adjust the length of the Abduction and stabilizing belts, you must move the adjusters on the front belt straps up or down.
- Use. In order to put the child on the Abduction and stabilizing belts, put the child on the installed belts, then wrap the child's leg and fasten both belt buckles to both front straps.

In order to remove the Abduction and stabilizing belts from the child, unfasten both belt buckles from the front straps.



Fig. 11.3.2



Fig. 11.3.3

#### 11.4. Seat minimizing kit

The Yeti stroller is equipped with two pairs of two inserts, two right and two left (Fig. 11.4). They are used to adjust the inside dimensions of the seat width. The upholstery of the seat minimizing kit is removable thanks to the zipper sewn in, which allows it to be replaced or cleaned.

- Assembly. The inserts are installed anywhere in the stroller, at any height adapting to the patient's indications. Installation is done with Velcro on the bottom, back and outer side of the insert.
- Disassembly is done by removing the element - pulling the velcro from the stroller upholstery.



Fig. 11.4

#### 11.5. The head restraint

The headrest is used to stabilize the user's head and is a softening element of the seat on the head section (Fig. 11.5). The headrest upholstery is removable thanks to the zipper sewn in, which allows it to be replaced or cleaned.

- Mounting the headrest with Velcro at the back of the headrest at the selected height.
- Disassembly is done by removing the element - pulling the Velcro from the stroller upholstery.



Fig. 11.5

#### 11.6. Adjustable canopy with window

The Yeti stroller can be equipped with an adjustable canopy with a window that helps protect the child from the sun or rain. The elements of the canopy are: canopy frame (front and rear arch), canopy opening adjustment mechanism (right and left), two snaps with red buttons fastening the canopy to the frame, canopy upholstery, and a mosquito net with a zipper.

- Assembly. The fastening of the adjustable canopy with the window are two plastic elements located on the outer sides of the upper part of the backrest. These attachments are used to mount the canopy with the window.



Fig. 11.6.1

To install the adjustable canopy with the window, place the backrest upright, and then slide both canopy latches into the backrest mounts in parallel. Mount the canopy in such a way that the upholstery at the back/bottom of the canopy overlaps the back of the stroller. In addition, make sure that the canopy adjustment element is secured with upholstery inside. This can

be done by stretching the upholstery and placing the loop on the button in the corner of the canopy. (Fig. 11.6.1)

- Disassembly. To disassemble the adjustable canopy with window, press the two red buttons on the roof fixing latches and slide the latches out of the backrest mounts.
- Use. To change the size of the canopy, hold the front arch and the rear arch of the canopy frame and bring the arches closer or closer to each other to set the desired position. The canopy surface adjustment mechanism latches automatically.

A special feature of the canopy is a window with a mosquito net for ventilation. (Rys.11.6.2). The window opens and closes with a zipper.



Fig. 11.6.2

### 11.7. Travel bag

The large and elegant travel bag is mounted on two handles on the inside of the handle to push the stroller (see Fig. 11.7)

The maximum bag load is 2kg.



Fig. 11.7

### 11.8. Shopping basket



Fig. 11.8

The Yeti stroller can be equipped with a shopping basket (Fig. 11.8.), which is used to transport small, ordinary items of everyday life, resistant to shocks, impacts, bruises, soaking, flooding and dirt.

It is mounted at the bottom of the stroller by wrapping it on the frame of the stroller and fastening it with a zipper.

The maximum load of the shopping basket is 3 kg.

### 11.9. Wedge

The Yeti stroller can be equipped with a wedge (Fig. 11.9.1.), which in the case of specific indications provides discharge (inclination to the sides) of the child's lower limbs from his body.

- **Assembly.** The wedge is mounted in the front of the seat, to the element connecting the upholstery with the barrier. (see Fig. 11.9.2). To do this, fasten the dedicated section of the upholstery belt on the barrier and put a wedge on it sticking Velcro together. This can be done on the front or back of the connecting strap, adjusting the distance to the size of the user. Then wrap the wedge and belt with a wrap, securing the wedge against slipping with Velcro.

- **Disassembly.** To remove the wedge, remove the wrapper and wedge and then the upholstery belt connecting it with the barrier detach and transfer to a specially prepared place in the seat upholstery.



Fig. 11.9.1



Fig. 11.9.2

### 11.10. Raincover

The Yeti stroller can be equipped with a rain cover that helps protect the child from moisture. This optional equipment requires mounting to other optional equipment. To mount the rain cover, put it on top of the entire adjustable canopy with window and then cover the seat and footrest with the cover (see Fig. 11.10).



Fig. 11.10

### 11.11. Winter gloves

Two winter gloves that increase the comfort of a child's assistant during colder weather can be put on the handle to push the stroller (see Fig.11.11). First you need to put on two material clamps and fasten them with velcro on the handle to push the stroller, then gloves are put on these material clamps and fastened with snaps.



Fig. 11.11



Fig. 11.12

### 11.12. Winter sleeping bag

The Yeti stroller can be equipped with a Winter sleeping bag, increasing the child's comfort during colder weather (see Fig. 11.12). This optional equipment requires installation to other optional equipment - a barrier.

To mount the Winter sleeping bag, place it on the footrest from below, and transfer it through the safety barrier from above and fasten it with snaps to the frame of the stroller.

### 11.13. Mosquito net

The Yeti stroller can be equipped with a mosquito net (see Fig. 11.13), which prevents insects from coming into direct contact with the child. This additional

accessory requires fitting to other optional equipment. In order to mount the mosquito net, put it from above on the entire adjustable canopy with window, and then cover the mosquito net seat and footrest.



Fig. 11.13

### 11.14. Sun umbrella

The Yeti stroller can be equipped with a Sun umbrella (see Fig. 11.14.1) that helps protect the child from the sun.

To install a sun umbrella, attach the umbrella holder to the trolley frame (see Fig. 11.14.2) and slide the Sun umbrella into the holder.



Fig. 11.14.1

Fig. 11.





## 12. Technical specification of the Yeti stroller

stroller length	120 cm
stroller width	66 cm
stroller height	134 cm
total basic mass	od 22kg
effective seat depth	31-45 cm
seat width	31*/43 cm
back height	65-85 cm
backrest width	30*/42 cm
distance between footrest and seat	29-44 cm
folded stroller length for transport	111 cm
folded stroller width for transport	66 cm
height of the folded stroller for transport	43 cm
basket load capacity	3 kg
bag capacity	2 kg
maximum load of the stroller	80 kg

The manufacturer reserves the right to make changes to technical parameters. Measurement values tolerance +/- 20mm, +/- 1.5kg, +/- 1°.

### 12.1. User parameters

max. patient weight	75kg
patient height	130-175 cm

### 12.2. Material composition of the stroller

Powder coated steel and aluminium, fabrics 100% polyester and plastic (PP, PCV, ABS, PA, PE, PU).

## 13. Delivery and tools

The Yeti stroller is delivered in a cardboard box and pre-assembled. The original packaging should contain the following items:

- Folded stroller frame with seat, backrest and wheels
- Optional equipment as ordered
- User manual with warranty card

Before use, check and make sure that the carton contains all the elements and none of them has been damaged, e.g. during transport. In the event that it is found that in the carton received is missing or damaged any component, contact the Seller and the Manufacturer in order to receive the missing elements or replace the product with a new, fully complete and functional.

The multifunctionality and extensive equipment of the Yeti stroller allow you to tailor the stroller to your personal requirements. To do this, read the entire manual carefully before use.

No tools are required to assemble and disassemble the Yeti stroller.

A Phillips screwdriver is required to adjust the height of the footrest. An Allen key is needed to adjust the backrest height and seat depth.

## 14. Assembly and disassembly of the stroller

- Stroller assembly. The product pre-packed in the original carton is already pre-assembled. Complete the installation by completing the following steps:
  - a) mount the rear wheels on the axle - see point 9.6
  - b) unfold the stroller frame with a footrest adjusting the parameters:
    - enable/disable lock of the direction of travel of the front wheels – see point 9.5
    - adjust the height and angle of the stroller push handle – see point 9.8
    - adjust the seat depth – see point 9.11
    - adjust the seat tilt angle – see point 9.12
    - adjust the height of the backrest – see point 9.13
    - set the tilt angle of the backrest – see point 9.14
    - adjust the length of the five-point belts – see point 10.4
    - adjust the height of the footrest and the tilt angle of the footrest – see point 9.15
  - c) adjust the upholstery – see points 10.1, 10.2, 10.3
  - d) install appropriate accessories and adjust parameters:
    - adjust the usable width of the seat (additional option) – see point 11.4
    - adjust the headrest position (additional option) – see point 11.5
    - install a safety barrier or a therapeutic table (additional options) – see points 11.1, 11.2
    - adjust the deflection and stabilization belts (additional option) – see point 11.3
    - adjust the canopy with the window (additional option) – see point 11.6
    - install and adjust the wedge (additional option) – see point 11.9
    - put on a travel bag (additional option)– see point 11.7
    - put on a shopping basket (additional option) – see point 11.8
    - put on selected items of equipment dedicated to the seasons:
      - Winter sleeping bag (additional option) – see point 11.12
      - winter gloves (additional option) – see point 11.11
      - mosquito net (additional option) – see point 11.13
      - raincover (additional option) – see point 11.10
      - umbrella (additional option) – see point 11.14
- Dismantling the stroller. In order to dismantle the stroller it is necessary to:
  - a) remove the optional equipment
  - b) unclip the upholstery (side elements) from the side frames of the stroller
  - c) lock the front wheels, fold the footrest and stroller frame
  - d) unclip the rear wheels
  - e) fold in the stroller closing it with the transport lock

## 15. Preparation for use

Before using the stroller:

- it is necessary to ensure that the child's footwear is comfortable and stabilizes the ankle well.
- make sure that the child is not placed in the stroller without clothing covering the body in places of contact with the upholstery of the stroller. This is to prevent chafing of the child's body and hygiene of the child.
- For children who have problems with urine retention, put a special oilcloth on the seat of the stroller, which can be purchased in pharmacies or medical stores.

- When using any additional aids such as venflons, drains, bags, catheters, etc., make sure that they are well secured against falling out, pulling out, clamping and do not interfere when using the stroller.

In order to start using the stroller properly, all the operations listed below must be carried out:

- Check the good technical condition of the stroller and all its components, and if the technical condition of the stroller or its individual component is bad or the user has any doubts about it, then the stroller must not be used. Please contact the Seller and the Manufacturer.
- Place the stroller on a flat and secure surface.
- Unfold the stroller correctly (see point 9.2)
- Engage the central brake.
- Make all necessary adjustments to the stroller (see chapter 9 and 10)
- Prepare for five-point belts, straps on the footrest and Abduction and stabilizing belts, if these accessories are available.
- Place the child in the stroller and make sure that all adjustments previously made are appropriate and safe for the child.
- Attach the child to five-point belts and straps on the footrest.
- Put on the child other necessary stroller equipment, if available.

## 16. Using the stroller

### 16.1. Carrying the stroller

The best way to carry a stroller is to use it with its wheels. The manufacturer allows the possibility of driving the folded stroller also with the use of wheels in the position indicated in Fig. 16.1.

In cases where this is not possible (e.g. to transport the stroller in a car), for safety reasons, the stroller can be lifted and moved after dismantling the additional equipment and the rear wheels, in order to reduce the load.

When carrying or lifting the stroller frame with mounted wheels, extreme caution should be exercised because the weight of this element exceeds 23kg. It is recommended that the stroller be carried by 2 people.

When carrying or lifting the stroller frame with the wheels mounted, fold the stroller frame and then hold the stroller frame only by the frame's fixed elements.

For safety reasons, do not lift the entire stroller, i.e. the frame with the seat and backrest, wheels, footrest and mounted accessories, as this may damage your health.

Do not lift the stroller by its footrest or safety barrier, therapy table or any other accessories, as this may damage the product and endanger the safety of others.

### 16.2. Getting on and off the stroller

When getting on and off the stroller, place the stroller as close as possible to a stable and safe place where the user will change. After making sure that the central brake is on and the footrest is raised to avoid standing on it, the user can change.



Fig. 16.1



### 16.3. Correct position in a stroller

For safety and health reasons, it is recommended that the user have the pelvis and back as close as possible to the backrest as much as possible to support the spine evenly (optimal support for all parts of the spine so that there is no free space).

The thighs should be supported on the largest possible surface, along their entire length, which is used to adjust the depth of the stroller seat (see point 9.11), and the feet should be supported on the footrest plate (see point 9.15), which prevents slipping, maintains optimal leg muscle tension and carries some of the weight from the upper body.

### 16.4. Riding in a stroller

When using the stroller, always hold the handle to push the stroller to prevent the wheel from rolling slowly. Release the central brake while holding the push handle of the stroller and push the stroller in a safe direction.

### 16.5. Moving on slopes

Always position the stroller forwards when climbing or driveways. Limit the speed on sloping surfaces and move as slowly as possible. The stroller assistant (parent, guardian) must have sufficient strength to control the stability of the stroller.

### 16.6. Overcoming thresholds or curbs

When overcoming small obstacles (small thresholds, small curbs, etc.), block the tilting foot with the foot and the handle to push the stroller downwards and towards you to raise the front wheels of the stroller and carefully drive the front wheels over the obstacle. Then slowly pull the push handle upwards away from the ground to raise the rear wheels of the stroller and carefully drive the rear wheels onto the obstacle.

### 16.7. Transport in a car and other vehicles

- To transport the Yeti stroller in the car, remove all accessories from the stroller and store them in a safe place. Then dismantle the stroller (see chapter 14) and place it in the trunk.

If the dismantled stroller is not in the trunk of the car separated from the passenger compartment, then the stroller frame must be securely attached to the vehicle. For this purpose, unused seat belts may be used.

- In order to transport a stroller in other vehicles, such as trains, trams, buses, planes, ships, etc., please ask the persons responsible for transport in these means of transport about the conditions of safe travel to the destination with a stroller and the appropriate requirements.
- In order to transport the stroller on other vehicles, remove all accessories from the stroller and store them in a safe place. Then dismantle the stroller (see chapter 14) and place it in the baggage storage area of the given vehicle. If the dismantled stroller is not in the trunk of the vehicle separated from the passenger compartment, the stroller frame must be securely attached to the vehicle.

The Yeti stroller should not be used as a seat to transport a child in a vehicle.

## 16.8. Cleaning and care

Small leaks and stains can be removed from the surface with a kitchen sponge and using clean, warm water.

Stains should be removed immediately after they appear on the surface of the stroller. Do not rub - rubbing will spread the stain to a larger surface and may distort the texture of the fabric.

The treatment of water stains such as stains from food, drinks and vegetable oils is as follows:

- Dry the stain by pressing a white paper towel, handkerchief or absorbent clean cloth over its entire surface.
- Soak a small piece of cloth with a weak solution of detergent and warm water (or cold water in the event of blood staining) and rinse the stain by removing it from the outside to the inside to prevent it from growing. Do not use a too wet cloth, as excessive amount of liquid may cause ripples after drying or enlargement of the stain and damage of the upholstery or foam inside.
- It may be necessary to repeat the above operation depending on the type of stain and the amount of stained material. The above operation should be continued until the stain is completely removed, but always remember to remove excess moisture from the stained material before rinsing again.
- Then rinse the material using only lukewarm water without detergent, using a clean sponge or a clean white cloth. As above, rinse the material in place after the stain in the right direction, i.e. from the outside in.

Disinfectants that do not contain organic solvents can be used in accordance with the manufacturer's instructions.

Do not dip the fabric. Make sure the product is always dry before use.

Do not use organic solvents or dry cleaning agents.

For disinfection, you can use a small amount of ordinary household bleach, diluted with water with a maximum concentration of 3% bleach.

ATTENTION! The manufacturer is not liable for damages resulting from the use of inappropriate cleaning products.

## 17. Storage of the stroller

The stroller should be stored in a clean, dry and well ventilated room, preferably at 0-35°C.

### 17.1. Preparation of the stroller for long-term storage

In order to prepare the Yeti stroller for long-term storage, secure the product and perform all service operations in the "Periodic Inspections" chapter as if they had occurred.

The stroller should be stored in a clean, dry and well ventilated room.

### 17.2. Corrosion protection

To avoid corrosive damage to the stroller, it should not be used or left in wet and damp rooms.

When the truck is wet, wipe all wet parts with a soft cloth and dry with a light oil preservative. This does not apply to the upholstery of the stroller, the provision of cleanliness and maintenance of which is described in detail in the Cleaning and care chapter.

As a standard, the stroller is covered with the highest quality coatings. In the event of damage to the varnish, protect the damage by applying a touch-up compound using the varnish proposed by the varnish Manufacturer.

Regular cleaning and lubrication of the stroller mechanisms with light oil guarantees long and trouble-free use of the stroller.

## 19. Re-use

The product is suitable for reuse.

All reusable products are subject to a particular burden.

Therefore, their features and performance may deviate from the new product to such an extent that the safety of patients and third parties during their use is at risk.

The product with proper care and maintenance retains its quality for a period of 3 years and this is the expected period of its use.

For secondary use, please:

- Remove all upholstered parts and replace them with new ones.
- The stroller frame with wheels and steering elements should be thoroughly cleaned and disinfected with commercially available cleaning and disinfecting agents.
- The product must then be checked by an Authorized Service Center (see Chapter 24) or the Manufacturer for wear and damage.
- All worn parts and components that do not match or are unsuitable for the user should be replaced with new ones.
- The product prepared in this way can be put back into use.

## 20. Service and periodic inspections

For safe and trouble-free use of the Yeti stroller, it must be regularly and thoroughly inspected and serviced as indicated in the table below.

In a situation where the stroller user is not able to carry out the activities himself, he should return the stroller to a specialized stroller service point or directly to the Manufacturer, who will properly review and will perform necessary service activities. These activities may be charged extra.

All repairs should be carried out in specialized, authorized Services (see point) or directly at the Manufacturer.

Only original spare parts that are supplied or recommended by the Manufacturer may be used in the stroller.

### Service should be performed:

- **Each time before using the stroller**
  - Check the efficiency of the central brake - After activating the central brake, both rear wheels of the stroller should be locked and cannot rotate.
  - Check the attachment of structural elements - It is necessary to check all mechanisms of the stroller and its movable elements, in particular, pay attention to: the correct mounting of the therapeutic seat with the backrest to the stroller frame, the correct mounting of all four wheels, the correct functioning of the five-point seat belts, properly locked seat angle adjustment mechanism, properly locked backrest angle adjustment mechanism and footrest correct mounting.
  - Check the condition of the tires - Check that the tires on all four wheels are not damaged or worn.
  - Inspection of solid frame parts, seat with backrest and footrest frame - Visually inspect the condition of the structural profiles. These elements must not be damaged or cracked.

- **Every one month**
  - Cleaning the stroller - Clean any moving parts, bearings, adjustment mechanisms with a soft cloth and then lubricate them with a few drops of light oil. Be careful not to soil the upholstery of the stroller or your own clothes.
  - Check all bolted connections for looseness - Carefully check that all bolts and knobs in the stroller are properly tightened.
- **Every six months**
  - Check the cleanliness of the upholstery and inspect the general condition of the stroller - Do not allow the stroller to be neglected as described in the chapter „Cleaning and care”.


## 21. Disposal

Disposal and recycling: the used stroller should be returned to a recycling center in accordance with the regulations in force in the given country.


## 22. Compliance with requirements for medical devices

The Yeti stroller meets the essential requirements of Council of Europe Directive 93/42 EEC amended by 2007/47/EC for medical devices. The Yeti stroller, in accordance with Directive 93/42 EEC, is a non-invasive, inactive medical device of class I according to rule 1.

- Declaration of Conformity:



MYWAM Kupiec, Bartold, Angres spółka jawna  
ul. Szczecińska 10, 41-516 Chorzów, Polska  
tel. 0048 32 733 11 31  
www.mywam.eu



---


**DEKLARACJA ZGODNOŚCI WE**  
**EC DECLARATION OF CONFORMITY**

<u>Wytwórca:</u>	<b>MYWAM Kupiec, Bartold, Angres spółka jawna</b>
<u>Manufacturer:</u>	<b>ul. Szczecińska 10, 41-516 Chorzów, Poland</b>
<u>Produkt:</u>	<b>Wózek inwalidzki specjalny YETI</b>
<u>Product:</u>	<b>Specialist Stroller YETI</b>

*Deklaruję z pełną odpowiedzialnością, że wózek inwalidzki specjalny YETI, będący wyrobem medycznym klasy I, jest zgodny z wymaganiami dyrektywy Rady Europy o wyrobach medycznych - 93/42/EEC z 14 czerwca 1993 roku zmienionej przez 2007/47/EEC z 21 marca 2010 roku, przyjętej w Polsce na mocy ustawy o wyrobach medycznych z dnia 20 maja 2010 r. (Dz. U. nr 107, poz. 679), oraz zgodny z normami zharmonizowanymi:*


*I solely declare that the specialist Stroller YETI is a class I device, and it is in conformity with the requirements of EU Council Directive for Medical Devices 93/42/EEC of 14 June 1993 amended by 2007/47/EEC of 21 March 2010, implemented in Poland under the Medical Devices Act of 20 May 2010 (Dz. U. nr 107, poz. 679), and is in conformity with harmonized standards:*

PN-EN 12183:2014	PN-EN 12182:2012	ISO 7176-1:2014	ISO 7176-3:2012
ISO 7176-5:2008	PN-ISO 7176-7:2001	ISO 7176-8:2014	PN-EN 1021-1:2007
PN EN 980:2010	PN EN 1041:2010	PN EN 14971:2011	




Chorzów, 14.10.2019 / 14 OCTOBER 2019

Data/Date



Podpis/Signature



**MYWAM**  
Kupiec, Bartold, Angres Spółka Jawna  
41-516 Chorzów, ul. Szczecińska 10  
www.mywam.pl  
tel. 32/733 11 31, mail: mywam@mywam.pl  
NIP: 627 273 67 93, Regon: 243428050  
KRS: 0000487471

### 23. Manufacturer

MyWam Kupiec, Bartold, Angres Spółka Jawna  
ul. Szczecińska 10, 41-516 Chorzów, Poland  
contact phone number: +48 32 733 11 31  
Website: [www.mywam.eu](http://www.mywam.eu)  
E-mail address: [mywam@mywam.pl](mailto:mywam@mywam.pl)




**The Yeti stroller is a Polish product. Polish quality, Polish price!**

### 24. Notes

### 25. List of authorized service points

POLAND

MyWam Kupiec, Bartold, Angres spółka jawna  
ul. Szczecińska 10, 41-516 Chorzów, Poland  
tel: +48 32 733 11 31  
Website: [www.mywam.pl](http://www.mywam.pl)  
E-mail adress: [mywam@mywam.pl](mailto:mywam@mywam.pl)



## 26. Warranty Card

Series and number: .....

Date of purchase: ..... - ..... - .....

Seller stamp and signature: .....

.....

## 27. Warranty conditions

- The manufacturer grants the buyer a warranty for the purchased stroller, for a period of 24 months from the date of purchase.
- Defects revealed during this period will be removed free of charge.
- The Buyer submitting a complaint on the stroller should deliver the stroller to the Seller or to the Manufacturer in the safety packaging, together with the elements of the stroller protected against damage.
- Damage caused by insufficient securing of the stroller shipment and its components is not subject to warranty repairs.
- The complainant may not demand gratification if he improperly marks or incorrectly insures the shipment of the stroller and its components.
- The cost of shipping the complained stroller to the Seller or to the Manufacturer shall be borne by the buyer complaining the stroller.
- The warranty period is extended by the time needed to repair the stroller.
- The activities provided for in the manuals, which the user is obliged to perform on their own and at their own expense, cannot be performed as part of a warranty repair.
- The buyer has the right to replace the complained stroller with a new stroller if the same damages or defects that have already been removed twice and which prevent the stroller from being used are still present.
- The stroller should be used in accordance with these operating manual. Unauthorized use of the stroller, inconsistent with these operating manual or inconsistent with the principles of correct operation will void the warranty.
- Making any structural changes to the stroller will void the warranty.
- Damage to the stroller caused by: improper operation, improper maintenance, incorrect storage of the stroller, mechanical damage as well as natural abrasion and wear of parts (e.g. upholstery, wheels) are not covered by the warranty. They are repaired at the user's expense.
- The manufacturer warns against performing warranty repairs at unauthorized service facilities. Finding such repair or unauthorized alteration of the stroller will void the warranty.
- The warranty card is valid provided that it has: the date of sale entered, confirmed by the Seller's company seal, the Seller's signature, type and serial number of the stroller.
- The warranty card is valid as long as it does not contain deletions and corrections.
- The stroller manufacturer runs its paid post-warranty service.
- In matters not covered by this warranty card, the provisions of the Civil Code shall apply.
- All matters related to the warranty should be directed to the Seller or directly to the Manufacturer.