

PRODUCER / ADVERTISER: MyWam LLC

CE



Illustrative photo



Operation and maintenance manual

Edition no.: 10 / 2023-12-05

The special stroller - **YETI**

The special pediatric stroller - **YETI YOUNG**

This is a medical device.

Use it according to the instruction manual or label.

The user with visual, reading or cognitive disabilities is requested to contact the Manufacturer.

Purpose of the Manual

The operation and maintenance manual is intended for those who provide direct care for children with disabilities (user assistants - parents, caregivers) who will be transported in the stroller. It is also intended for doctors and physical therapists operating the stroller.

ATTENTION: Before using the stroller, each User is required to read and follow the contents of this manual.

NOTE: In the event of any medical incident related to the Yeti Special Stroller and Yeti Young Special Stroller, the Manufacturer should be notified immediately, in writing, by contacting via registered letter or email with acknowledgement of receipt.

Manufacturer contact information:

- postal address: MyWam LLC, Lwowska 34 Street, 41-500 Chorzów, Poland
- contact phone number: +48 32 733 11 31
- e-mail address: mywam@mywam.pl



Version of the Operation and Maintenance Manual: EN **Edition 10 / 2023-12-05**

MyWam LLC reserves the right to make technical and commercial changes to the contents of the Operation and Maintenance Manual without prior notice. Each change will be identified by the date of the last update of the manual. **The** electronic version of the manual will be available at www.mywam.eu

NOTE: The pictures contained in the User Manual are for reference only and may differ from the purchased product.

Table of contents

1. Purpose and indications for use of Yeti/ Yeti Young stroller.....	5
2. Quick start	5
▪ Security.....	6
3. Important safety instructions and rules	6
4. Identification plate	7
5. Warning lables.....	7
▪ Product	8
6. General construction of the stroller	8
7. Configuration and use – stroller frame	9
7.1. Frame folding mechanism.....	9
7.2. Frame unfolding mechanism.....	9
7.3. Folding the stroller	10
7.4. Transport lock. Engaging and disengaging of the transport lock.....	10
7.5. Front swivel wheels. Front wheels direction lock, mounting, dismounting.....	10
7.6. Rear wheels. Mounting and dismounting of the rear wheel on the axle	11
7.7. Rear wheel suspension system.....	11
7.8. Push handle. Mechanism and angular adjustment.	11
7.9. Central brake. Engaging and disengaging of the central brake	11
7.10. Tilting foot.....	11
7.11. Seat depth adjustment.....	11
7.12. Adjusting the seat angle (tilt).	12
7.13. Backrest height adjustment.....	12
7.14. Adjusting the backrest recline angle.....	12
7.15. Footplate. Adjustment and use.	12
8. Configuration and use - therapeutic seat with backrest	13
8.1. Stroller upholstery. Installation of upholstery to the frame of the stroller frame. Removal.....	13
8.2. Adjustable tension system	14
8.3. Seat and backrest upholstery adjustment.....	14
8.4. Five-point belts.....	15
9. Additional equipment.....	15
9.1. Safety barrier.....	15
9.2. Therapy table	15
9.3. Abduction and stabilizing belts.....	15
9.4. Seat minimizing kit.....	16
9.5. Headrest	16

9.6. Adjustable canopy with window	16
9.7. Travel bag.....	17
9.8. Shopping basket.....	17
9.9. Wedge (abductoin block).....	17
9.10 Raincover.....	17
9.11 Winter gloves	18
9.12 Winter sleeping bag	18
9.13 Mosquito net.....	18
9.14. Sun umbrella	18
10. Technical specifications of the Yeti stroller	19
10.1. Material composition of the stroller	19
11. Delivery and tools	19
12. Assembly and disassembly of the stroller	20
13. Preparation for use.....	21
14. Using the stroller	21
14.1. Carrying the stroller	21
14.2. Getting on and off the stroller.....	21
14.3. Correct position in a stroller.....	22
14.4. Riding in a stroller.....	22
14.5. Moving on slopes.....	22
14.6. Overcoming thresholds or curbs	22
14.7. Transport in a car and other vehicles	22
14.8. Cleaning and care	23
15. Storage of the stroller.....	23
16. Service and periodic inspections	23
17. Disposal.....	24
18. Compliance with requirements for medical devices	24
List of authorized service points.....	Błąd! Nie zdefiniowano zakładki.
20. Warranty Card	25
21. Notes.....	26

1. Purpose and indications for use of Yeti/ Yeti Young stroller

The Yeti Special Stroller and Yeti Young Special Baby Stroller is a Class I medical device indicated for use by those disabled persons who, due to their condition, cannot independently walk stand and/or sit down. The product is intended for the transportation and movement of the above-mentioned persons in a sitting or lying position, as well as for postural correction, and this is its intended use.

Clinical indications include cerebral palsy, meningeal hernia, dystrophies, paresis, etc.

The product is intended for users with a height of 130-175cm and a weight of up to 80kg. The product should be used as recommended by the attending physician or physiotherapist.

Contraindications:

In case of a contraindication to sitting and/or passive verticalization, the adjustment of the stroller should be consulted with a physiotherapist or attending physician.

2. Quick start

In order to use the stroller quickly, it should be unpacked and prepared for use according to paragraph 12 Mounting and dismounting the stroller and paragraph 13 Preparing for use.



▪ Security

3. Important safety instructions and rules

- It is forbidden to use the stroller in a manner for which the stroller is not intended.
- The product contains small parts that can be dangerous if swallowed by the user!
- Due to the high adjustability of the stroller, during the angular adjustment of the seat and the angular adjustment of the backrest, it is possible to set a position in which the user's head is below the pelvic line. It is forbidden to use such a position, as it could lead to the danger of the user swallowing or choking.
- It is forbidden to place additional weight on the stroller, which especially applies to placing additional weight on the back of the backrest, such as a backpack or bag. This may cause the stroller to tip over and injure the user! Does not apply to original additional accessories dedicated to the Yeti stroller not exceeding the allowable load. All additional items attached to the stroller affect its stability.
- **NOTE:** Depending on the adjustments made - refers to the setting of the stroller's push handle, backrest height adjustment, backrest angle adjustment, seat depth adjustment, seat angle adjustment, footrest angle adjustment and footrest setting - this stroller may exceed the recommended maximum dimensions.

- **CAUTION:** When passing through doors or similar constrictions that restrict the direction of travel of the stroller, make sure there is enough space on both sides of the stroller to prevent hands or other body parts from getting pinched, as well as to prevent damage to the stroller.
- **NOTE:** The temperature of the stroller may increase when the stroller is exposed to sunlight.


WARNING!



- Never leave the user unattended in the stroller, as there is a risk of injury and there is a risk of accident!
 - This product is not suitable for running or skating or rollerblading.
 - The stroller is not intended as a seat in motor vehicles.
 - Beware of gaps and moving parts - risk of jamming.
 - Do not move the stroller with the user.
- The stroller is packaged in a bag. Keep out of reach of infants and small children.

- **CAUTION:** It is forbidden to use the stroller, if the technical condition of the product or its single component raises doubts in the user's mind or is clearly abnormal, then it must not be used. In such a situation, it is recommended to urgently contact the Seller for advice.


4. Identification plate



UD code plate



Warning plate



Nameplate

5. Warning labels

<p>Product serial number</p>	<p>Manufacturer name</p>
<p>Permitted user total weight</p>	<p>The product is a medical device</p>
<p>Production date</p>	<p>DESCRIPTION "WARNING!" and the symbol in this manual are used to enhance the reader's attention to the content marked with this sign. Non-compliance with content marked with this description may endanger the health and safety of the user and accompanying persons, and may also damage the product.</p>
<p>Conformity assessment mark according to the applicable European regulations and directives</p>	
<p>Warning about proceeding in accordance with the product Operation and Maintenance Manual</p>	

▪ **Product**

6. General construction of the stroller

• **Steel stroller frame including seat and backrest**



- 1- Handle for pushing the stroller - see section 7.8
- 2- String adjustment of upholstery excess - see paragraph 8.3
- 3- Five-point belts - see paragraph 8.4
- 4- Removable back and seat upholstery - see point 8.1
- 5- Slot for retracting excess upholstery - see paragraph 8.1
- 6- Frame wrapper - see section 8.1
- 7- Adjustable footrest with straps - see paragraph 7.15
- 8- Tilting foot - see paragraph 7.10
- 9- Central brake - see section 7.9
- 10- Rear wheel suspension system - see section 7.7
- 11- Rear wheels, solid - see paragraph 7.6

- 12- Angle adjustment mechanism of the push handle-see paragraph 7.8
- 13- Frame folding mechanism - see section 7.1
- 14- Backrest angle adjustment - see paragraph 7.14
- 15- Seat depth adjustment - see paragraph 7.11
- 16- Seat angle adjustment - see paragraph 7.12
- 17- Backrest height adjustment - see point 7.13
- 18- Locking mechanism of the frame folding mechanism - see section 7.3
- 19- Angle adjustment mechanism of the whole footrest - see paragraph 7.15
- 20- Backrest angle adjustment mechanism - see point 7.14
- 21- Transport lock - see section 7.4
- 22- Footrest length adjustment - see paragraph 7.15
- 23- Rear wheel quick release mechanism - see section 7.6
- 24- Locking the direction of travel of the front wheels - see paragraph 7.5
- 25- Front swivel wheels, solid - see paragraph 7.5

• **Accessory equipment**

- Safety barrier-see paragraph 9.1
- Therapy table-see paragraph 9.2
- Abduction and stabilization belts-see section 9.3
- Tapering inserts-see section 9.4
- Adjustable headrest ⇕- see paragraph 9.5
- Adjustable canopy with window-see paragraph 9.6
- Travel bag-see paragraph 9.7
- Shopping basket-see paragraph 9.8

- Wedge (abductoin block)-see paragraph 9.9
- Raincover-see paragraph 9.10
- Winter gloves-see paragraph 9.11
- Sleeping bag for legs-see paragraph 9.12
- Mosquito net-see paragraph 9.13
- Sunshade umbrella-see paragraph 9.14

7. Configuration and use – stroller frame

7.1. Frame folding mechanism

- right and left breaker/joint - a component that allows the upper frame to be folded to reduce its dimensions,
- Frame folding mechanism including locking mechanism.

7.2. Frame unfolding mechanism

Lay the frame of the stroller horizontally on the ground and pull down the transport lock (Fig. 7.2.1). Then gently pull the backrest to an upright position



Fig. 7.1

(Fig. 7.2.2).



Fig. 7.2.1 Fig



Fig. 7.2.2

To set up the stroller, pull the handle by lifting it gently upwards while placing the stroller in a position close to the desired one (Fig. 7.2.3) until the frame folding lock is engaged (Fig. 7.2.4). The frame of the stroller will be correctly unfolded only when the frame folding mechanism is properly closed. This is evidenced by a distinct "click" sound.



Fig. 7.2.3



Fig. 7.2.4

7.3. Folding the stroller

• **Fold the frame:** release the lock of the frame folding mechanism (Fig. 7.3) by pulling it upward, and pull both frame folding handles upward until the frame is released. Fold the frame of the stroller by folding it forward, towards the front wheels. Make sure the push handle and footrest are in a position that does not interfere with each other. Pull the backrest release mechanism. Put on the transport lock.



Fig. 7.3

7.4. Transport lock. Engaging and disengaging of the transport lock

- **Engaging the lock:** fold the frame, put the lock eyelet on the lock pin so that two overlap tightly (Fig. 7.4).
- **Disengaging the lock:** before unfolding the frame, pull the lock eyelet off the lock pin and only then can you proceed to unfold the frame of the stroller.



Fig. 7.4

the



Fig. 7.5.1

7.5. Front swivel wheels. Front wheels direction lock, mounting, dismounting

- **Directional lock:** insert the lock knob into the slot (Fig. 7.5.1).
- **Pulling off the lock:** pull off and turn the lock knob and leave it out of the slot (Fig. 7.5.2).



Fig. 7.5.2

7.6. Rear wheels. Mounting and dismounting of the rear wheel on the axle



Fig. 7.6.1

Mounting the wheels on the axle: press the quick release button (Fig. 7.6.1) and put the wheel on the beginning of its axle. Then stop pressing on the quick release button and slide the wheels onto its axle until it clicks.

Removal of wheels from the axle: press the quick release button located in the center of the wheel and slide the out of its axle (Fig. 7.6.2).



Fig. 7.6.2

wheel

7.7. Rear wheel suspension system

The shock absorber operates automatically. It is adjusted by activating the updown lever located on the bottom of the shock absorber (Fig. 7.7). When the lever is switched, the shock absorber spring will be compressed (hard ground option) or stretched (soft ground).




Fig. 7.7



Fig. 7.8

7.8. Push handle. Mechanism and angular adjustment.

- **Adjustment mechanism:** two buttons on the outside of the frame, located on the joints of the push handle of the stroller (Fig. 7.8).
- **Angle adjustment:** press both buttons of the angle adjustment mechanism of the push handle at the same time and change its angle of tilt, then release them.

7.9. Central brake. Engaging and disengaging of the central brake

- **Engaging:** press the brake pedal down (Fig. 7.9), towards the ground, until it clicks. The center brake is used only for total wheel locking - there is no dynamic braking function!
- **Disengaging:** pull the brake pedal up with your foot.



Fig. 7.9

7.10. Tilting foot

Using a tilting foot makes it easier to drive the stroller over small thresholds or curbs (Fig. 7.10).



Fig. 7.10

7.11. Seat depth adjustment



Attach the central brake, loosen the bolt with an allen wrench on both sides of the frame (Fig. 7.11) and slide out the seat frame (seat frame) evenly into the desired position. Secure the correct position by tightening the screw on the sides of the frame. The frame can be slid out leaving a min 15cm overlap of the steel frame components. The frame is adjustable from 30 to 45 cm. The adjustment range of the seat upholstery is described in section 8.3.

7.12. Adjusting the seat angle (tilt)

- **Seat angle adjustment mechanism:** engage the central brake, unscrew the knob from both sides, adjust the seat according to the selected position and screw the knob in the correct position (Fig. 7.12).



Fig. 7.12



Fig. 7.13

7.13. Backrest height adjustment

- **Adjustment:** engage the central brake, loosen the screw with an allen wrench on both sides of the frame (Fig. 7.13), and slide or slide the backrest frame (backrest frame) evenly to the desired position.

Secure the selected position by tightening the screw. The frame can be extended leaving a min 15cm overlap of the steel frame elements. The maximum level is indicated by a ball lock. The frame is adjustable from 65 to 85 cm.

7.14. Adjusting the backrest recline angle

- **Adjustment:** engage the central brake, pull the lock release cable upwards (Fig. 7.14.2) and adjust the backrest to the desired value (8 degree scale).



Fig. 7.14.1

Lock the backrest by releasing the lock release cable and after making sure that the lock has snapped on both sides of the stroller at the same time (went deep into the ridges).

When changing the position of the backrest, it is necessary to adjust the upholstery.

The backrest angle adjustment mechanism changes the position of the backrest relative to the seat therapeutic in six positions, ranging from a sitting position at right angles to an almost fully reclined position.




Fig. 7.14.2

Note: Every time you change the angle of the back rest, you need to readjust the length of the five-point straps.

7.15. Footplate. Adjustment and use.

- **Adjust the height of the footrest:** unscrew the knobs securing the footrest to the frame (Fig. 7.15.1), on both sides of the footrest, adjust the footrest to the desired height and screw it back in with both screws, tightening them with the knobs until you feel resistance.

- **Angle adjustment mechanism for the entire footrest:** two plastic black joints with buttons in the center (Fig. 7.15.2). Press the buttons on both sides, set the desired position and release the buttons.



Fig. 7.15.1

The mechanism for angular adjustment of the footrest during use should be protected by upholstery wraps (see paragraph 8.1).



Fig. 7.15.2

- **Use of the footrest:** it is recommended to set the footrest plate at such a height that the user's feet can freely rest on it with their entire surface, while at the same time the user's thighs are parallel to the plane of the seat. In the footrest plate there are straps securing the user's feet. It is recommended to put them on when using the stroller.

By default, the footrest plate is in the horizontal position. However, when transporting, storing, rehabilitating, getting in and out of the stroller, the footrest plate must be folded and rotated towards the footrest frame.

8. Configuration and use - therapeutic seat with backrest

8.1. Stroller upholstery. Installation of upholstery to the frame of the stroller frame. Removal.

Upholstery consists of several components: upholstery base, seat upholstery, backrest upholstery and frame wraps. Before adjusting the upholstery, adjust the frame accordingly: backrest height, seat depth and footrest height.

- **Upholstery assembly:**

- ✓ Fix the upholstery **base** in place of the seat with adjustable tapes (Velcro) stretching them carefully. Install the tapes of the upholstery base from the bottom by wrapping them on the seat frame (from the outside). Install the first Velcro on the footrest side still on the plastic part of the breaker/joint, but from the inside of the seat. Wrap the excess upholstery in the minimum setting of the seat on the back of the seat. The base of the upholstery should be properly stretched allowing, moreover, to put the wrappers on the frame at the seat breakers/joints (at the contact with the footrest, on the right and left sides).
- ✓ **Frame wraps** protect the user from friction at breakaway points. Mount with Velcro downward on plastic breakaway pieces.
- ✓ Place **the** upholstery of the **backrest** on the upholstery base and install with the straps from the back wrapping the frame from the outside. If the minimum height of the backrest is used, wind the excess upholstery from the top onto the frame and pull it down with the straps adjusting it to the correct distance through the lower buckle. In the case where the backrest is extended in the maximum the upper part put on the frame and pull down through the upper buckle.

To properly install the side upholstery covers of the backrest of the stroller, it is recommended to position the backrest to the 45° position and fasten with Velcro on the right and left side of the stroller in two places on the side frame (at the folding lock of the frame). Once the backrest is in the upright position, the excess upholstery from the sides can be collected with cords with clamps.

- ✓ Apply **the seat upholstery** to the end by properly arranging and adjusting the Velcro on the upholstery base. Pull the excess upholstery under the backrest upholstery and wrap upward (over the back) under the backrest upholstery straps. The seat upholstery mounting strap (mounted under the bottom of the seat covering the frame from the outside) should fit over the first Velcro of the base upholstery and should not interfere with the frame wrapper.

Velcro the bottom of the seat upholstery to the footrest frame.



The seat upholstery of the stroller has an additional strip of material connecting it to the railing. In case it is not needed, it can be hidden by pulling it in the dedicated slot (Fig. 8.1).

Fig. 8.1

- **Upholstery removal:** dismantle using the previous steps in reverse order, so:
 - ✓ Unfasten the seat upholstery straps and pull this upholstery off by "detaching" it from the Velcro straps,
 - ✓ Unhook the sides of the backrest in the backrest upholstery and pull off the protection of the upper part of the backrest,

- ✓ Unhook all the straps and pull off the backrest upholstery,
- ✓ Pull the wraps off the frame,
- ✓ Open all the upholstery base straps and pull them off the frame.

8.2. Adjustable tension system



Fig. 8.2.1

The adjustable tension system is a set of as many as 14 independent velcro straps that allow the tension of both the seat and the backrest to be adjusted to suit the individual needs of the user. This system consists of 7 straps on the back of the backrest (Fig. 8.2.1) and 7 straps on the bottom of the seat (Fig. 8.2.2).

- **Adjustable:** tighten or loosen the straps as appropriate. Possibility of point or all-over collapse in the seat and backrest, and increase support after by all-over or point increase in stiffness/tension of both seat and backrest.

8.3. Seat and backrest upholstery adjustment

- **Adjust seat upholstery:** slide out the frame (see paragraph 7.11) and adjust the upholstery to the required parameters. To change the seat depth, remove the seat



Fig. 8.2.2

upholstery, adjust the seat frame for depth (slide the seat frame in or out to the desired size). Put the seat upholstery on the seat base and fasten it with Velcro.

- **Backrest upholstery adjustment:** slide out the frame (see paragraph 7.13) and adjust the upholstery to the required parameters. To change the height of the backrest, unhook the upper part of the backrest upholstery from the clamping straps on the back of the stroller and adjust the upper frame of the stroller (slide out or slide

in the upper part of the backrest frame until it reaches the correct position). After selecting the position, arrange the upper part of the upholstery adjusting it to the

frame. With the stroller set in the minimum backrest position, you may need to wind the upholstery onto the top of the frame. Fasten the upholstery from the top with straps to the buckles (lower buckles for the stroller with the smallest backrest, upper buckles for the seat in the largest position). Before inserting the user, adjust the tension of the seat and backrest upholstery by adjusting the set of straps (See section 8.2).



Fig. 8.3.1 Fig. 8.3.2

- To arrange the sides of the backrest upholstery, adjust the upholstery's "pull-down" state with the cords. For a stroller in the sitting position, pull the excess upholstery with strings in two places on each of the sides (Fig. 8.3.1). On the other hand, before adjusting the backrest to the reclining position, the cords should be completely loosened releasing the possibility of maximum unfolding of the upholstery (Fig. 8.3.2).

8.4. Five-point belts

The belts have been factory-fitted to the stroller and should not be removed, as they are used to carry the user safely, prevent the user from slipping, and are also helpful for additional stabilization of the user's pelvis and torso (Fig. 8.4.1).

- **Adjustment:** slide the adjusters on the harness up or down and move the protectors. To adjust the height of the strap attachment to the backrest, move the harness attachment from the low attachment to the high attachment or vice versa.
- **Use:** put the harness on the user's right and left sides and push both harnesses into the belt buckle (Fig. 8.4.2), which is fixed on the seat. To tighten the belts press the button on the top (chest) buckle and the bottom buckle, release both from the buckle and pull the harnesses off the user's right and left sides. **Note:** Every time you change the angle of the back rest, you need to readjust the length of the five-point straps.



Fig. 8.4.1 harnesses



Fig. 8.4.2

9. Additional equipment

9.1. Safety barrier

The safety barrier is mounted with the same fasteners as the therapy table, so you can't use both at the same time.

- **Installation:** two plastic pieces (snaps) mounted at the tapered point of the frame (Fig. 9.1). Open the two clamps located on the railing and simultaneously slide them onto the frame and snap them into place. The railing will be properly installed only when the railing goes around the frame snapping into place. This is evidenced by a distinct "click" sound.

- **Disassembly:** pull off both railing latches located on the frame and at the same time slide out the railing.



Fig. 9.1

9.2. Therapy table

The maximum load of the table is 3kg.

- **Installation:** two plastic pieces (snaps) mounted at the constriction of the frame (Fig. 9.2). Open the two clamps (latches), simultaneously slide them onto the frame and snap them into place. The table will be properly installed only when the latch goes around the frame snapping into place.

This is evidenced by a distinct "click" sound.

- **Disassembly:** pull off both latches located on the frame and at the same time slide out the table.



Fig. 9.2

9.3. Abduction and stabilizing belts



Fig. 9.3.1

Abduction and stabilization straps (Fig. 9.3.1) are used for additional stabilization of the pelvis and to prevent the user from slipping out of the stroller. Fasten on the user's thighs so that the user's pelvis is as close to the back rest as possible.

- **Installation:** in four places - attach the two rear straps to the rear frame of the seat with buckles (Fig. 9.3.2), attach the two front straps to the frame of the stroller under the narrowing of the railing (sideways) putting the belt through the factory loop in the strap (Fig. 9.3.3)



Fig. 9.3.2



Fig. 9.3.3

- **Adjustment:** slide the adjusters on the front belt straps up or down.
- **Use:** To put on the belts seat the user on the installed belts, wrap around the user's leg and fasten both belt buckles to both front straps. To remove the belts, undo both belt buckles from the front straps.

9.4. Seat minimizing kit

The inserts are used to adjust the internal dimensions of the seat width. The upholstery of the inserts is removable thanks to a sewn-in zipper, making it possible to replace or clean it.

- **Mounting:** we mount anywhere in the stroller, at any height adapting to the patient's indications. Mounted with Velcro on the bottom, back and outside of the insert.
- **Disassembly:** peel off the Velcro from the upholstery of the stroller.



Fig. 9.4

9.5. Headrest

The headrest serves to stabilize the user's head and is an element that softens the seat at the head section (Fig. 9.5). The upholstery of the headrest is removable thanks to a sewn-in zipper, making it possible to replace or clean it.

- **Installation:** mount with Velcro on the back of the headrest at the desired height.
- **Disassembly:** peel off the Velcro from the upholstery of the stroller.



Fig. 9.5

9.6. Adjustable canopy with window

- **Canopy components:** canopy frame (front and rear arch), canopy opening adjustment mechanism (right and left), two snaps with red buttons securing the canopy to the frame, canopy upholstery, and mosquito net window with zipper.

- **Installation:** the mounts are two plastic pieces located on the outer sides of the top of the back rest. Set the backrest in an upright position, slide the two canopy snaps to the mounts on the backrest.



Fig. 9.6.1

Make sure that the upholstery inside the canopy secures the canopy adjustment element (pull the upholstery and put the loop on the button in the corner of the canopy - Fig. 9.6.1).

- **Disassembly:** press the two red buttons on the canopy mounting tabs and slide the tabs out of the backrest mounts onto the back.
- **Use:** To change the size of the canopy surface, grasp the front arch and the back arch of the canopy frame, and by bringing the two arches closer or further apart or toward each other, adjust the desired position. The canopy surface size adjustment mechanism will snap into place automatically.

The canopy has a window with a mosquito net for ventilation. (Fig. 9.6.2). The window opens and closes with a zipper.



9.7. Travel bag



The bag is mounted on two handles, located on the inside of the handle to pushing the stroller (Fig. 9.7). The maximum load of the bag is 2kg.

Fig. 9.7

9.8. Shopping basket

The basket is mounted at the bottom of the stroller by wrapping it around the frame of the stroller and fastening it with a zipper.

The maximum load of the shopping basket is 3kg.



Fig. 9.8

9.9. Wedge (abductoin block)

The wedge (abductoin block), for certain indications, provides an abduction (side-to-side deviation) of the user's lower extremities from their body center (Fig. 9.9.1).



Fig. 9.9.1

- **Installation:** fasten a section of the upholstery belt on the railing and to it apply a wedge (abductoin block) by gluing the Velcro together. You can install at the front or back of the connecting belt, adjusting the distance to the size of the user. Then wrap the wedge (abductoin block) and belt with Velcro securing the wedge (abductoin block) from slipping (Fig. 9.9.2).

- **Disassembly:** pull off the wrap, wedge (abductoin block) and upholstery belt connecting it to the railing, detach it and put it in a specially prepared place in the seat upholstery.



Fig.9.9.2

9.10 Raincover

Requires mounting to another accessory - a canopy.

Installation: put the film from the top on the canopy and cover the seat of the stroller the footrest with the film (Fig. 9.10).

9.11 Winter gloves



Fig. 9.10

- **Assembly:** put on two fabric clamps, fasten them with Velcro on the push handle of the stroller, put on the gloves and fasten (Fig. 9.11).

Requires mounting to another accessory - railing (Fig. 9.12).

- **Assembly:** put the sleeping bag from the bottom on the footrest, and from the top put it through the safety railing and clip it to the frame of the stroller.



Fig. 9.11

9.12 Winter sleeping bag

Requires mounting to another accessory - railing (Fig. 9.12).



Installation: put the sleeping bag from the bottom on the footrest, and from the top put it through the safety railing and clip it to the frame of the stroller.

9.13 Mosquito net

Requires mounting to another accessory - a canopy (Fig. 9.13).

- **Installation:** apply the mosquito net from the top of the canopy and drape the mosquito net over the stroller seat and footrest.



Fig. 9.13

Fig. 9.12.

9.14. Sun umbrella



Fig. 9.14.1

- **Installation:** mount the umbrella holder to the stroller frame (Fig. 9.14.2) and slide the umbrella into this holder.



Fig. 9.14.2

10. Technical specifications of the Yeti stroller

length of the stroller	114.5 - 146 cm
the width of the stroller	66cm
the height of the stroller	117-138 cm
total weight	24 - 28.6 kg
effective seat depth	31-45 cm
seat width	28* / 43 cm
height of the backrest	65 - 85cm
backrest width	25*/42 cm
distance of the footrest from the seat / number of adjustments	28.5 - 43.5 cm / 7 positions
length of the stroller, folded up for transport	110 cm
width of the stroller, folded up for transport	66cm
height of the stroller, folded up for transport	50 cm
carrying capacity of the basket	3 kg
carrying capacity of a travel bag	2 kg
max. patient weight	80 kg
patient height	130-175 cm
seat plane angle / number of adjustments	approx. 9.5°-29° / 5
backrest angle / number of adjustments	approx. 3°-99° / 8 + lying down
footrest angle / number of adjustments	about 80°-180° /6 positions

NOTES:

*) values given with inserts

Tolerance of measurement values +/- 20mm, +/- 1.5kg, +/- 1°

10.1. Material composition of the stroller

Powder-coated steel and aluminum, 100% polyester fabrics and plastic (PP, PVC, ABS, PA, PE, PU).

11. Delivery and tools

The original packaging should include the following items:

- Folded stroller frame with seat and backrest and wheels,
- Accessory equipment according to the order,
- User manual with warranty card,
- 10 mm Allen wrench,

Before use:

- a. check for damage to the packaging (if the packaging is damaged, unpack in the presence of the supplier and check that none of the components have been damaged),
- b. check in accordance with the specifications the contents of the shipment (if any of the components are missing, contact the Seller).

No tools are required to assemble and disassemble the stroller.

A Phillips screwdriver is needed to adjust the height of the footrest. An allen wrench is needed to adjust the height of the backrest and the depth of the seat.

12. Assembly and disassembly of the stroller

- Stroller assembly. The product factory-packed in the original stroller on is already pre-assembled.

Complete the assembly by performing the following steps:

- a) Install rear wheels on the axle - see fig. 7.6
 - b) unfold the frame of the stroller together with the footrest by adjusting the parameters:
 - switch on / off the directional lock of the front wheels - see paragraph 7.5
 - adjust the height and angle of the push handle of the stroller - see section 7.8
 - adjust the seat depth - see paragraph 7.11
 - adjust the seat tilt angle - see section 7.12
 - adjust the height of the backrest - see paragraph 7.13
 - adjust the backrest recline angle - see paragraph 7.14
 - adjust the length of the five-point belts - see section 8.4
 - adjust the height of the footrest and the angle of the footrest - see paragraph 7.15
 - c) adjust upholstery - see 8.1, 8.2, 8.3
 - d) Install the appropriate accessories and adjust the parameters:
 - adjust the usable width of the seat (additional option:) - see section 9.4
 - adjust the position of the headrest (additional option) - see section 9.5
 - install a safety barrier or a therapy table (additional options) - see 9.1, 9.2
 - adjust the abduction and stabilization straps (additional option) - see section 9.3
 - adjust the canopy with window (additional option) - see paragraph 9.6
 - Fit and adjust the wedge (abductoin block) (optional extra) - see section 9.9
 - put on the travel bag (additional option) - see paragraph 9.7
 - put on a shopping basket (additional option) - see paragraph 9.8
 - put on selected pieces of equipment dedicated to the seasons:
 - sleeping bag for legs (additional option) - see paragraph 9.12
 - winter gloves (additional option) - see paragraph 9.11 - mosquito net (additional option) - see paragraph 9.13
 - Raincover (additional option) - see paragraph 9.10
 - umbrella (additional option) - see paragraph 9.14
- Stroller disassembly. In order to disassemble the stroller, it is necessary to:
 - a) dismantle the accessory
 - b) detach the upholstery (side elements) from the side frames of the stroller
 - c) lock the front wheels, fold the footrest and the frame of the stroller

- d) detach the rear wheels
- e) Fold the stroller by locking it with a transport lock.

13. Preparation for use

Before using the stroller, you should:

- Check that the frame of the stroller is properly unfolded and that the frame folding mechanism is properly closed (locked).
- Check that the seat has been properly mounted to the frame of the stroller
 - Correctly make all adjustments and fit the stroller to the user's needs, including adjusting the proper depth of the therapy seat and the height of the backrest, the proper tilt of the therapy seat and the backrest recline, adjust the proper angle of the footrest, adjust the proper tension of the straps on the backrest and on the seat, check the length of the five-point straps. Check that all adjustments are properly seated and secured.
- Check that the upholstery of the stroller is dry.
- Secure the user with all available stabilizing elements, such as five-point belts, narrowing inserts, abduction and stabilizing belts, wedge (abductor block), footrest, footrest straps, safety barrier or therapy table. This is to prevent the user from falling and slipping out of the stroller and to prevent the user's feet from falling off the footrest. Make sure that all adjustments made beforehand are appropriate and safe for the user.
- Ensure that the user is not placed in the stroller without clothing to cover the body in areas of contact with the upholstery of the stroller. This is to prevent abrasions to the user's body and for hygienic reasons for the user.
- For users with continence problems, put a special oilcloth on the seat of the stroller, which can be purchased at pharmacies or medical stores.
- If you use any additional aids such as venflons, drains, bags, catheters, etc., make sure that they are well protected against falling out, pulling out, clamping, and do not interfere with the use of the stroller.

14. Using the stroller

14.1. Carrying the stroller

The best way to move a stroller is to drive it using its wheels (Fig. 14.1).

In cases where this is not possible (for example, to transport the stroller in a car), then for safety reasons, the stroller can be lifted and carried with both hands after removing the accessories to reduce the weight to be carried and folding them. Special care should be taken when moving or lifting the frame of the stroller. The frame should be moved after it is folded and the transport lock is in place. Do not use accessory equipment to move the stroller - you risk damaging

it.



Fig.14.1

14.2. Getting on and off the stroller

When getting in and out of the stroller, place the stroller as close as possible to a stable and safe place for the user to transfer to. After making sure the rear brakes are engaged and the footrest has been removed to avoid standing on it, the user can be placed in the stroller.

Do not stand on the footrest plate, as there is a risk of tipping over the stroller, injury to the user and risk of damage to the product.

Make sure that when getting in or out of the stroller, the child and no other person should stand on the footrest, as this risks damaging it and injuring the user!

14.3. **Correct position in a stroller**

For safety and health reasons, it is recommended that the user have the pelvis and back positioned as close to the backrest as possible, as far as possible, in order to support the spine evenly (optimal support of all parts of the spine, so that no free space is created).

The thighs should be supported as much as possible, along their entire length, which is helped by adjusting the depth of the seat of the stroller (see paragraph 7.11), and the feet should be rested on the footrest plate (see paragraph 7.15), which prevents sliding, maintains optimal muscle tension in the legs and transfers some of the weight from the higher parts of the body. **NOTE:** Due to the high adjustability of the stroller, during the angular adjustment of the seat and the angular adjustment of the backrest, it is possible to set such a position in which the user's head is below the pelvic line. It is forbidden to use such a position, as it could lead to the danger of the user swallowing or choking.

14.4. **Riding in a stroller**

When using the stroller, always hold the stroller's push handle to prevent the stroller from rolling inertly. Release the center brake while holding the push handle of the stroller and push the stroller in a safe direction.

14.5. **Moving on slopes**

Always position the stroller facing forward when climbing hills or ramps. Limit your speed on inclined surfaces and move as slowly as possible. The person driving the stroller must have sufficient strength to control the stability of the stroller.

Always keep in mind that the position of the center of gravity of the stroller changes when traversing uneven or sloping terrain. In such a case, be very careful and belay the stroller.

When using driveways, make sure that the driveway is suitable, that the wheels of the stroller are high enough to negotiate the driveway, that the frame of the stroller will not come into contact with the driveway, and that a second person is available to assist, as there is a risk of the stroller tipping over, damaging the product and injuring the user.

14.6. **Overcoming thresholds or curbs**

When overcoming small obstacles (small thresholds, small curbs, etc.), press down the tilt foot and the push handle of the stroller downward and toward you with your foot to lift the front wheels of the stroller and carefully drive the front wheels over the obstacle. Then slowly pull the push handle of the stroller upward from the ground to lift the rear wheels of the stroller and carefully drive the rear wheels onto the obstacle.

14.7. **Transport in a car and other vehicles**

- Remove all accessories from the stroller and store them in a safe place. Fold the stroller and place it in the trunk.

If the dismantled stroller is not in the trunk of the car separate from the passenger compartment, the stroller frame must be carefully secured to the vehicle. For this purpose, you can use unused seat belts.

- In order to transport a stroller in other vehicles, such as trains, streetcars, buses, airplanes, ships, etc., it is necessary to ask the people in charge of transportation in these means of transport about the conditions for safe travel to the destination with a stroller and the relevant requirements.



Yeti family strollers should not be used as a seat to carry a child in a vehicle.

14.8. Cleaning and care

Upholstery can only be hand washed at up to 30°.

Stains should be removed as soon as they form on the surface of the stroller.

CAUTION: Do not clean with chemicals containing solvents and bleach. The manufacturer is not responsible for damages caused by the use of inappropriate cleaning products.

Make sure the product is always dry before use.

When using detergents dedicated to this type of material - follow the instructions of the detergent used.



Hand wash only



Do not bleach



Do not tumble dry
drum



Dry in position
horizontal



Do not iron

15. Storage of the stroller

Clean stroller should be stored in a clean, dry and well-ventilated room, preferably at 0-35 degrees. For long-term storage, additionally follow the steps described in point. 16 "every one month".

16. Service and periodic inspections

For safe and trouble-free use of the stroller, it is necessary to regularly and thoroughly subject it to the inspections and servicing described below.

In a situation where the stroller user is not able to carry out the above-mentioned activities on their own, they should return the stroller to a specialized stroller service center or directly to the Manufacturer, who will properly carry out the inspection and necessary servicing activities. These activities may be subject to an additional fee.

NOTE: The Manufacturer is responsible for defects or technical damage to the stroller caused only when the stroller is operated correctly. Responsibility for the proper operation of the stroller rests with the Manufacturer only if all recommendations of the Manufacturer are strictly followed.

NOTE: The manufacturer is responsible for the proper functioning of the stroller only if original parts and original accessories are used in it.

All repairs must be carried out at specialized authorized Service Centers (see paragraph 19) or directly at the Manufacturer. Repairs to the stroller performed by persons other than an authorized Service Center or the Manufacturer will void the warranty. The Manufacturer is not responsible for any injury or damage to the product that may be a consequence of repairs performed by persons other than authorized Service or the Manufacturer.

Maintenance activities should be performed:

- **Each time before using the stroller**

- Checking the efficiency of the central brake - When the central brake is applied, both rear wheels of the stroller should be locked and must not rotate.
- Checking the attachment of structural components - Inspect all mechanisms of the stroller and its moving parts, and in particular, pay attention to: proper installation of the therapeutic seat with backrest to the

stroller frame, proper installation of all four wheels, proper functioning of the fivepoint safety belts, properly locked seat angle adjustment mechanism, properly locked backrest angle adjustment mechanism, and proper installation of the footrest.

- Checking the condition of the tires - Check the tires of all four wheels for damage or wear.
- Inspection of the fixed components of the frame, seat with backrest and footrest frame - Visually inspect the condition of the structural profiles. These components must not be damaged or cracked.
- **Every one month**
 - Cleaning the stroller
 - Visual inspection of moving parts, bearings, adjustment mechanisms - lubricate with silicone grease.
 - Check all bolted connections for play.

17. Disposal

Disposal and recycling: in accordance with country-specific regulations.

18. Compliance with requirements for medical devices

Class I Medical Device in accordance with the Regulation of the European Parliament and of the Council (EU) 2017/745 of April 5, 2017 on medical devices. The manufacturer has issued an appropriate declaration of conformity for the product, which can be found at www.mywam.eu

19. List of authorized service points

POLAND
MyWam
Wiosny Ludów 91 street, 40-373 Katowice, Poland
tel: +48 32 733 11 31
www.mywam.eu
adres e-mail: serwis@mywam.pl



LLC



20. Warranty Card

Series and number:.....

Purchase date: - -

Seller's stamp and signature:

.....
.....
.....

Warranty terms and conditions

- The seller grants the buyer a warranty for the purchased stroller, for a period of 24 months, starting from the date of purchase.
- **The buyer complaining about the stroller should deliver the stroller to the Seller in its packaging, clean (free of stains, additional equipment not belonging to the stroller, moisture, hair, etc.) along with the stroller components protected from damage.**
- **Damage caused by insufficient protection of the shipment of the stroller and its components is not subject to warranty repairs.**
- The claimant may not claim gratification if he mislabeled or improperly insured the shipment of the stroller and its components.
- The warranty period is extended by the time the stroller is repaired.
- The activities provided for in the instructions, which the user is required to perform on his own and at his own expense, cannot be performed under warranty repair.
- Use of the stroller not in accordance with the intended use, not in accordance with this manual or not in accordance with the rules of proper operation will void the warranty.
- Making any structural changes to the stroller will void the warranty (in the particular area of change or in its entirety).
- Damage to the stroller caused by: improper operation, improper maintenance, improper storage of the stroller, mechanical damage and natural wear and tear of parts (e.g. upholstery, wheels) are not covered by the warranty.
- The manufacturer warns against making warranty repairs at unauthorized service facilities. The discovery of such repair or arbitrary alteration of the stroller will void the warranty.
- **The warranty card is valid provided that it has: the inscribed date of sale confirmed by the Seller's company stamp, the Seller's signature, the type and factory number of the stroller.**
- The warranty card is valid provided that it does not contain deletions and corrections.
- The manufacturer of the forklift provides post-warranty service at its headquarters for a fee.
- In matters not covered by this warranty card, the provisions of the Civil Code (Poland) shall apply.
- All matters related to the warranty and complaints of the stroller, please contact the Seller.

21. Notes