



POLISH MANUFACTURER
OF SPECIAL NEEDS STROLLERS

C € MD MEDICAL DEVICE - CLASS I

The MEWA special needs stroller is designed for children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for children at a height of 90-125 cm and a weight of 40 kg.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust MEWA to the required position and ensure the user with optimal protection and stability.

The shock-absorbing suspension on all four, large, gel-filled wheels together with the tilt-foot facilitate travel on uneven terrain. Modern and stylish design is an additional advantage of this stroller.

#### INDICATIONS

paralysis and paresis of lower limbs of various origins, paralysis and paresis of the trunk, weakening of postural muscles responsible for maintaining the correct body position, cerebral palsy, myelomeningocele, muscular dystrophies.

#### CONTRAINDICATIONS:

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

### colour options.



UPHOLSTERY



**FRAME** 

spray painted / glossy



UPHOLSTERY



**FRAME** spray painted / half matt













**FRAME** spray painted / glossy



**UPHOLSTERY**synthetic



**FRAME** spray painted / glossy



UPHOLSTERY synthetic



FRAME spray painted / glossy













**NOTE**: The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist.

**NOTE**: Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.

### features. safety and comfort

**TILT IN SPACE SEAT** 

simple mechanism



5,5°-21,5°

4 positions

**RECLINING BACKREST** wide range of adjustments



8°-86°

4 positions

### ADJUSTABLE AND ANGLE TILT FOOTREST

simple and quick adjustment



64°- 180°

5 positions

#### **SEAT MINIMIZING KIT**

reducing inserts



The SEAT MINIMIZING KIT in MEWA is a set of 4 comfortable, independent inserts that allows you to reduce both width and depth of the seat.

#### STABILIZATION SYSTEM

adjustable pads



The head and torso support consists of two pairs of independent pads that are height, width and angle adjustable. Its wide range of adjustments and simplicity gives you a fast and effective adaptation to the user's individual needs.

#### LYING POSITION

wide range of backrest recline



With a fully adjustable seat and footrest you may quickly alter the stroller into a cozy bed perfect for a short rest or a nap.

#### **REVERSIBLE SEAT**

arrange your travel easily



The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time

#### **FOLDING LIGHTWEIGHT FRAME**

ease of storage, transportation and dismantling



The construction of main elements allows you to dismantle MEWA into 3 separate units (frame, seat and wheels).

#### **DRIVING SUPPORT**

comfort for the bumpiest journeys



The stroller is equipped with a shock-absorbing suspension on all four wheels and an adjustable rear suspension in the rear wheels, swivel front castor wheels with option to lock, central brake and a tilt-foot.

# functional accessories. in practice



## accessories.wide range of equipment



### parameters.

#### TECHNICAL DATA OF A STROLLER READY TO DRIVE



seat depth 24\*/28 cm

ubackrest height 62 cm

footrest adjustable within 17 - 32 cm

- total weight 20 28 kg
- total lenght 94 126 cm
- total height 124 129 cm
- total width 65 cm

### dimension after folding: frame L 84 x W 64 x H 26 cm / seat L 88 x W 43 x H 42 cm

#### **USER PARAMETERS**

- 🙀 growth 90 125 cm
- maximum weight 40 kg