

Pegaz

SPECIAL NEEDS STROLLER FOR CHILDREN

GENERAL INFORMATION

The **PEGAZ special needs stroller** is designed for the children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for users at a height of **80 – 140 cm** and a maximum weight of **40 kg**.

PEGAZ is an ultralight umbrella type stroller made of aluminium and steel alloys equipped with the **"I Grow With You"** system, narrowing inserts, reclining backrest and a wide range of accessories. Simple to fold and unhold, conveniently compact (starting from 10 kg) and easy to manoeuvre – **PEGAZ** will be your day-to-day support for many years.



  **MEDICAL DEVICE - CLASS I**

COLORS



gray upholstery



red upholstery



THIS IS MEDICAL DEVICE.
USE IT IN ACCORDANCE WITH THE INSTRUCTIONS FOR USE OR THE LABEL.

FEATURES / SAFETY AND COMFORT

SEAT POSITION angle parameter



25°

1 positions

BACKREST RECLINE angles parameters



22°, 34°

2 positions


ADJUSTABLE FOOTREST high regulation



20, 30, 40cm

3 positions

I GROW WITH YOU extensive adjustment options



FOLDING LIGHTWEIGHT FRAME ease of storage and transportation



The one-piece, **umbrella construction** of **Pegaz** stroller allows for quick and easy folding which due to its compact size and low weight facilitates daily use, storage and transport.

DRIVING SUPPORT comfort and maneuverability



The stroller is equipped with **shock-absorbing suspension** on all four wheels, **swivel front castor wheels** with option to lock, **two independent rear wheel brakes** and a **tilt-foot**.

INDICATIONS

- paralysis and paresis of lower limbs various origins
- paralysis and paresis of the trunk
- weakness of postural muscles, responsible for maintenance correct body position
- cerebral palsy
- myelomeningocele
- muscular dystrophies

CONTRAINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

FUNCTIONAL ACCESSORIES IN PRACTICE



EQUIPMENT



**NARROWING
INSERTS**



HEADREST



ABDUCTION BELTS



SAFETY BARRIER



FIVE-POINTS BELTS



**ADJUSTABLE
ABDUCTION BLOCK**



THERAPY TRAY



**CANOPY WITH
WINDOW**



SHOPPING BASKET



TRAVEL BAG



SUN UMBRELLA



RAIN COVER



SLEEPING BAG



MOSQUITO NET

TECHNICAL PARAMETERS OF THE STROLLER

seat width: 30*/38 cm

effective seat depth: 25/33 cm

backrest height: 50/60 cm

backrest width: 30*/37 cm

distance between footplate and seat: 20, 30, 40 cm, 3 positions

*values are given with inserts

total weight: 10 – 19,6 kg

stroller length: 90 – 108 cm

stroller height: 106 – 130 cm

stroller width: 56 cm

IMPORTANT NOTE! The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist. Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.