

Grizzly

GENERAL INFORMATION

The **GRIZZLY** special needs stroller is designed for the children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for children at a height of **85 – 140 cm** and a maximum weight of **45 kg**

GRIZZLY is equipped with the **"I Grow With You"** system that allows you for a smooth extension of the backrest height and seat depth as you child grows. Moreover **GRIZZLY** offers you a comprehensive stabilization system consisting of adjustable pads at the level of the head, hips and chest – together with a fully adjustable footrest, tilt in space seat and reclining backrest we bring you a stroller that is tailored to fit your child's needs and will serve for years.

Crash-tested for safety, Grizzly is perfectly suited for use in specially adapted vehicles.



  **MEDICAL DEVICE - CLASS I**

INDICATIONS

Paralysis and paresis of lower limbs of various origins, paralysis and paresis of the trunk, weakening of postural muscles responsible for maintaining the correct body position, cerebral palsy, myelomeningocele, muscular dystrophies.

CONTRINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.



THIS IS MEDICAL DEVICE.

USE IT IN ACCORDANCE WITH THE INSTRUCTIONS FOR USE OR THE LABEL.

FEATURES / SAFETY AND COMFORT

TILT IN SPACE SEAT

simple mechanism



0° – 30°

4 positions

RECLINING BACKREST

wide range of adjustments



91° – 212°

6 positions

ADJUSTABLE AND ANGLE TILT FOOTREST

simple and quick adjustment



55° – 225°

5 positions

I GROW WITH YOU

reducing insert



„I GROW WITH YOU” system allows to extend **the seat depth up to 37 cm**, **backrest height up to 70 cm** (no tools required).
Designed to adapt and evolve!

STABILIZATION SYSTEM

adjustable pads



A complete support system consisting of a **four independent pairs of pads (head, torso, hips and lateral)** guarantees maximum of adjustment possibilities as well as highest comfort and safety during travel.

LYING POSITION

wide range of backrest recline



With a fully adjustable seat and footrest you may quickly alter the stroller into cozy bed, perfect for a short rest or a nap.

REVERSIBLE SEAT

arrange your travel easily

FOLDING LIGHTWEIGHT FRAMEease of storage, transportation
and dismantling**DRIVING SUPPORT**

comfort for the bumpiest journeys



The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time.



The construction of main elements allows you to dismantle **Grizzly** into 2 separate units – the frame and the seat. All wheels are also removable.



The stroller is equipped with **an adjustable rear suspension** in the rear wheels, swivel front castor wheels with option to lock, **central brake** and **a tilt-foot**.

FUNCTIONAL ACCESSORIES IN PRACTICE

EQUIPMENT



SIDE PADS

HEAD/TORSO/HIPS
ADJUSTABLE PADS

ABDUCTION BELTS



SAFETY BARRIER



FIVE-POINTS BELTS

ADJUSTABLE
ABDUCTION BLOCK

THERAPY TRAY

CANOPY WITH
WINDOW

WINTER GLOVES



TRAVEL BAG



SHOPPING BASKET



SUN UMBRELLA



RAIN COVER



SLEEPING BAG



MOSQUITO NET



GRIZZLY BEAR SET

TECHNICAL PARAMETERS OF THE STROLLER

seat width: 18* – 40 cm
 effective seat depth: 25 – 37 cm
 backrest height: 57 – 70 cm
 backrest width: 18* – 34 cm
 distance between footplate and seat: 23 – 40,5 cm, 7 positions
 *values are given with pads

total weight: 20 – 24 kg
 stroller length: 82 – 135 cm
 stroller height: 108 – 123 cm
 stroller width: 65 cm

IMPORTANT NOTE! The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist. Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.