

# Pegaz

## SPECIAL NEEDS STROLLER FOR CHILDREN

### GENERAL INFORMATION

The **PEGAZ special needs stroller** is designed for the children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for users at a height of **80 – 140 cm** and a maximum weight of **40 kg**.

**PEGAZ** is an ultralight umbrella type stroller made of aluminium and steel alloys equipped with the **"I Grow With You"** system, narrowing inserts, reclining backrest and a wide range of accessories. Simple to fold and unhold, conveniently compact (starting from 10 kg) and easy to manoeuvre – **PEGAZ** will be your day-to-day support for many years.



  **MEDICAL DEVICE - CLASS I**

### COLORS



gray upholstery



red upholstery



**THIS IS MEDICAL DEVICE.**  
**USE IT IN ACCORDANCE WITH THE INSTRUCTIONS FOR USE OR THE LABEL.**

# FEATURES / SAFETY AND COMFORT

## SEAT POSITION angle parameter



25°

1 positions

## BACKREST RECLINE angles parameters



22°, 34°

2 positions

## ADJUSTABLE FOOTREST high regulation



20, 30, 40cm

3 positions

## I GROW WITH YOU extensive adjustment options



## FOLDING LIGHTWEIGHT FRAME ease of storage and transportation



The one-piece, **umbrella construction** of **Pegaz** stroller allows for quick and easy folding which due to its compact size and low weight facilitates daily use, storage and transport.

## DRIVING SUPPORT comfort and maneuverability



The stroller is equipped with **shock-absorbing suspension** on all four wheels, **swivel front castor wheels** with option to lock, **two independent rear wheel brakes** and a **tilt-foot**.



## INDICATIONS

- paralysis and paresis of lower limbs various origins
- paralysis and paresis of the trunk
- weakness of postural muscles, responsible for maintenance correct body position
- cerebral palsy
- myelomeningocele
- muscular dystrophies

## CONTRAINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

## FUNCTIONAL ACCESSORIES IN PRACTICE



# EQUIPMENT



**NARROWING  
INSERTS**



**HEADREST**



**ABDUCTION BELTS**



**SAFETY BARRIER**



**FIVE-POINTS BELTS**



**ADJUSTABLE  
ABDUCTION BLOCK**



**THERAPY TRAY**



**CANOPY WITH  
WINDOW**



**SHOPPING BASKET**



**TRAVEL BAG**



**SUN UMBRELLA**



**RAIN COVER**



**SLEEPING BAG**



**MOSQUITO NET**

## TECHNICAL PARAMETERS OF THE STROLLER

seat width: 30\*/38 cm

effective seat depth: 25/33 cm

backrest height: 50/60 cm

backrest width: 30\*/37 cm

distance between footplate and seat: 20, 30, 40 cm, 3 positions

\*values are given with inserts

total weight: 10 – 19,6 kg

stroller length: 90 – 108 cm

stroller height: 106 – 130 cm

stroller width: 56 cm

**IMPORTANT NOTE!** The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist. Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.