




CMD MEDICAL DEVICE - CLASS I

**NEW  
GENERATION**

# mewa

## SPECIAL NEEDS STROLLER FOR CHILDREN



## GENERAL INFORMATION

The **MEWA special needs stroller** is designed for the youngest users who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for children at a height of **90 – 125 cm** and a maximum weight of **40 kg**.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust **MEWA** to the required position and ensure the user with optimal protection and stability.

The Mewa stroller is equipped with rear-wheel suspension, which significantly improves comfort when navigating uneven terrain. This feature helps absorb shocks, ensuring a smoother ride for the child.

## INDICATIONS

- paralysis and paresis of lower limbs various origins
- paralysis and paresis of the trunk
- weakness of postural muscles, responsible for maintenance correct body position
- cerebral palsy
- myelomeningocele
- muscular dystrophies

## CONTRAINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

**THIS IS MEDICAL DEVICE.**

**USE IT IN ACCORDANCE WITH THE INSTRUCTIONS FOR USE OR THE LABEL.**



### LIGHTWEIGHT CONSTRUCTION

The lightweight design of the stroller is an ideal solution for parents who value mobility and comfort. The stroller is easy to carry and maneuver, which makes walks and trips more enjoyable. At the same time, it also makes it easier to fold and store the stroller.



### NEW UPHOLSTERY MATERIAL

The material used in the construction of the stroller is friendly to both people and animals. Made of safe, non-toxic raw materials, it meets the highest quality and ecological standards. It is gentle on the child's skin, and at the same time durable and resistant to damage.



### SHOCK ABSORPTION

Shock absorption in strollers is crucial for comfort and safety. It effectively dampens shocks on uneven terrain, which allows the child to travel in peace and without any discomfort. At the same time, it makes it easier to drive the stroller, allowing smooth maneuvering, which reduces the burden on the parent.



### BREATHABLE MATERIAL

The interior of the stroller and part of accessories are made of 100% breathable SONIC material that provides excellent air circulation. Thanks to this, the child feels comfortable even during long walks, and the skin remains dry and free from irritation. The material is not only breathable, but also soft and gentle, further enhancing user comfort.



### SHOPPING BASKET

The shopping basket with a capacity of up to 3 kg, mounted on the **Mewa** stroller, is a convenient solution for parents. It allows you to easily transport necessary accessories and small items, without the need to carry additional bags.



### DURABLE POLYURETHANE TIRES

Extremely durable and flexible tires, made entirely of polyurethane, so there is no need to pump them up. This design significantly reduces the risk of accidental damage, ensuring reliability and long service life.



### REVERSIBLE SEAT

The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time.



**NEW CANOPY DESIGN**

one-piece, with ventilation, easier to use

Redesign canopy blends elegance with functionality and offers superior sun protection, adapting to the weather. Built-in ventilation provides airflow, keeping your child cool on hot days.

**NEW BAG DESIGN****NEW UPHOLSTERY MATERIAL**

more durable and abrasion resistant

**UNIVERSAL ACCESSORIES****OTHER COLORS****BEIGE****GREEN****BLACK****GOLD**

# EQUIPMENT


**SET MINIMIZING KIT**

**HEAD/TORSO  
ADJUSTABLE PADS**

**ABDUCTION BELTS**

**SAFETY BARRIER**

**FIVE-POINTS BELTS**

**ADJUSTABLE  
ABDUCTION BLOCK**


**THERAPY TRAY**

**CANOPY WITH  
WINDOW**

**WINTER GLOVES**

**TRAVEL BAG**

**SHOPPING BASKET**

**SUN UMBRELLA**

**CUSHION PILLOWS**

**SLEEPING BAG**

**MOSQUITO NET**

**RAIN COVER**

## TECHNICAL PARAMETERS OF THE STROLLER

seat width: 25\*/32 cm  
 effective seat depth: 24 – 28 cm  
 backrest height: 62 cm  
 backrest width: 20\*/30 cm  
 distance between footplate and seat: 17 – 32 cm, 5 positions  
 \*values are given with pads

total weight: 20 – 28 kg  
 stroller length: 94 – 126 cm  
 stroller height: 124 – 129 cm  
 stroller width: 62,5 cm

**IMPORTANT NOTE!** The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist. Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.