

REDESIGNED FOR YOUNG USERS, BUILT TO ADAPT.

FOR USERS: 90 - 125 CM, UP TO 40 KG

NEW
GENERATION

mewa

The Mewa special needs stroller

Is designed for the youngest users who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust Mewa to the required position and ensure the user with optimal protection and stability.

The Mewa stroller is equipped with rear-wheel suspension, which significantly improves comfort when navigating uneven terrain. This feature helps absorb shocks, ensuring a smoother ride for the child.

CE **MD** MEDICAL DEVICE - CLASS I

**THIS IS A MEDICAL DEVICE.
USE IT ACCORDING TO THE INSTRUCTIONS FOR USE OR LABEL.**

PRODUCT SPECIFICATIONS

Lightweight construction

The lightweight design of the stroller is an ideal solution for parents who value mobility and comfort. The stroller is easy to carry and maneuver, which makes walks and trips more enjoyable. At the same time, it also makes it easier to fold and store the stroller.

Mesh material

Inside the stroller and in selected accessories, a mesh-structured material is used to support natural air circulation. This ensures that the interior of the stroller remains properly ventilated, and the user stays thermally comfortable even during longer use.

Reversible seat

The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time.



Regulations



SEAT ANGLE
ADJUSTMENT



BACKREST ANGLE
ADJUSTMENT



FOOTREST
ADJUSTMENT

6,5° - 22,5°

4 positions

9° - 87°

4 positions

64° - 180°

5 positions

TECHNICAL PARAMETERS



FOR USERS: 90 - 125 CM, UP TO 40 KG

Total weight: **20,5 - 25 kg**

Stroller length: **97 - 132 cm**

Stroller width: **62,5 cm**

Stroller height: **125 - 128 cm**

Effective seat depth: **24* / 28 cm**

Seat width: **25* / 32 cm**

Backrest height: **62 cm**

Backrest width: **20* / 30 cm**

Distance between a footrest and seat: **17 - 32 cm**

/ 5 positions

* values are given with inserts

FUNCTIONS



Adjustable Features

The Mewa stroller can be easily adapted to the child's height thanks to the seat minimizing kit, adjustable pads, abduction wedge, and adjustable footrest. As the child grows, these elements can be gradually removed or adjusted to maintain proper support and user comfort.



Rear wheels with shock-absorbing suspension

The rear wheels are equipped with a suspension system that effectively reduces vibrations when driving over uneven surfaces. This ensures greater comfort for the child and smoother stroller handling.



Swivel castor lock

The front wheels feature a directional lock, allowing for stable stroller guidance on straight paths and in more demanding terrain.



GREY

BEIGE

GREEN

BLACK / GOLD

COLOURS

ACCESSORIES



CANOPY WITH WINDOW



ADJUSTABLE PADS



INSERTS AND CUSHIONS



ABDUCTION BELTS



FIVE-POINTS BELTS



ABDUCTION WEDGE



TRAVEL BAG



SAFETY BARRIER



THERAPY TRAY



SHOPPING BASKET



WINTER GLOVES



SUN UMBRELLA



SLEEPING BAG



RAIN COVER



MOSQUITO NET