

# SAFE & COMFORTABLE MOBILITY FOR THE YOUNGEST.

FOR USERS: 80 - 110 CM, UP TO 30 KG

mouse

## The Mouse special needs stroller

Is designed for the youngest users who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust Mouse to the required position and ensure the user with optimal protection and stability.

The Mouse stroller features rear-wheel suspension and off-road wheels, making it much easier to travel over uneven terrain.

CE  MEDICAL DEVICE - CLASS I

**THIS IS A MEDICAL DEVICE.  
USE IT ACCORDING TO THE INSTRUCTIONS FOR USE OR LABEL.**

# PRODUCT SPECIFICATIONS

## Folding lightweight frame

The construction of the Mouse special needs stroller consists of two independent elements: the frame and the seat. The seat is attached to the frame using dedicated mounting adapters.



## Lying position

Thanks to the fully adjustable seat and footrest, the stroller can be easily converted into a lying position, providing comfortable support for rest or sleep. The extended, adjustable upholstery allows additional coverage, helping to protect the child from wind and external conditions.



## Reversible seat

The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time.



## Regulations



SEAT ANGLE ADJUSTMENT

8° - 23°

4 positions



BACKREST ANGLE ADJUSTMENT

16° - 88°

4 positions



FOOTREST ADJUSTMENT

84° - 173°

5 positions

# TECHNICAL PARAMETERS

# FUNCTIONS



## FOR USERS: 80 - 110 CM, UP TO 30 KG

Total weight: **17,2 - 23,2 kg**

Stroller length: **95 - 124 cm**

Stroller width: **60 cm**

Stroller height: **122 - 134 cm**

Effective seat depth: **22\* / 26 cm**

Seat width: **24\* / 32 cm**

Backrest height: **52 cm**

Backrest width: **30 cm**

Distance between a footrest and seat: **15 - 23 cm**

**/ 5 positions**

\* values are given with inserts

### Adjustable Features

The Mouse stroller can be easily adapted to the child's height thanks to the seat minimizing kit, adjustable pads, abduction wedge, and adjustable footrest. As the child grows, these elements can be gradually removed or adjusted to maintain proper support and user comfort.



### Rear wheels with shock-absorbing suspension

The rear wheels are equipped with a suspension system that effectively reduces vibrations when driving over uneven surfaces. This ensures greater comfort for the child and smoother stroller handling.



### Swivel castor lock

The front wheels feature a directional lock, allowing for stable stroller guidance on straight paths and in more demanding terrain.



### Redesigned footrest

This upgrade provides a wider range of height adjustment – now it can be set from 15 to 23 cm from the seat.

# BASIC EQUIPMENT



CANOPY WITH WINDOW



ADJUSTABLE PADS



INSERTS AND CUSHIONS



ABDUCTION BELTS



FIVE-POINTS BELTS



SAFETY BARRIER



TRAVEL BAG



SHOPPING BASKET



THERAPY TRAY



SUN UMBRELLA



WINTER GLOVES



ABDUCTION WEDGE



SLEEPING BAG



RAIN COVER



MOSQUITO NET

# ADDITIONAL ACCESSORIES